
INSTRUCTIONS FOR USE



Actim[®] Pancreatitis

EN	Instructions for use	3
DE	Gebrauchsanweisung	5
ES	Instrucciones de uso	7
FI	Käyttöohje	9
FR	Notice d'utilisation	11
IT	Istruzioni d'uso	13
NL	Gebruiksaanwijzing	15
PT	Instruções de utilização	17
RU	Инструкция по применению	19
SV	Bruksanvisning	21

Structure of dipstick

Test procedure and results

Instructions for use

Numbers ① - ⑧ refer to illustrations on inner cover.

Structure of dipstick

① Dip area ② Result area ③ Test line ④ Control line

Intended use

The Actim® Pancreatitis test is a visually interpreted, qualitative immunochromatographic dipstick test for use when acute pancreatitis is suspected. The test is based on detection of elevated levels of trypsinogen-2 in urine.

Kit components

Actim Pancreatitis kit 32731ETAC contains 10 and the kit 32732ETAC contains 20 dipsticks in sealed aluminum foil pouches with desiccant. Each kit contains also instructions for use.

Storage

Store the dipsticks at +2...+25 °C. Use the dipsticks shortly after their removal from the aluminum foil pouch.

Specimen collection

The specimen is urine. Urine samples should be tested as soon as possible after collection. Samples may be stored in a refrigerator at +2...+8 °C for 24 hours. If the sample needs to be stored for longer, it should be frozen. The sample should be brought to room temperature before testing.

Test procedure and interpretation of the results

1. If stored refrigerated, allow the aluminum foil pouch to reach room temperature. Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper pink part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch.
2. Place the yellow dip area into the sample ⑤ and hold it **there until you see the liquid front enter the result area ⑥**. Remove the dipstick from the sample and place it horizontal.
3. The result can be interpreted as positive as soon as two blue lines become visible in the result area. Negative result should be read at 5 minutes ⑦.
Do not pay attention to any lines appearing later than 5 minutes.

4. If **two blue lines, the test line and the control line**, appear, the test result is **positive**. If **one blue line, control line**, appears, the test result is **negative**.

If **the control line does not appear**, the test is **invalid ⑧**.

Limitations of the test

The test is intended for professional *in vitro* diagnostic use only.

Notes

- No quantitative interpretation should be made based on the test results.
- Care must be taken when placing the dipstick in sample. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing.
- Do not use the dipstick if its aluminium foil pouch is not intact.
- When dipping, be careful to hold the dipstick in position (with the dip area in the urine sample) until the sample liquid front reaches the result area
- If the control line does not appear, the test is invalid, and should be repeated using another dipstick.
- Excessive amount of biotin in urine sample may interfere for the formation of control line.
- If the test result cannot be interpreted clearly it is recommended that the test be repeated.
- At 5 minutes, the appearance of any faint-to dark blue test line along with a control line indicates a positive result. However, do not pay attention to any lines appearing after 5 minutes.
- If only the control line is visible, the result should be interpreted as negative only after 5 minutes have elapsed.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

Principle of the test

Levels of trypsinogen-2 have been shown to increase in the urine of patients with acute pancreatitis. The test is based on immunochromatography. It involves two monoclonal antibodies to human trypsinogen-2. One is bound to blue latex particles (the detecting label). The other is immobilized on a carrier membrane to catch labeled particles and indicate a positive result.

When the dip area of the dipstick is placed in a urine sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains trypsinogen-2 it binds to the antibody attached to the latex particles. The particles are carried by the liquid flow and, if trypsinogen-2 is bound to them, they bind to the catching antibody. A blue line (test line) will appear in the result area if the concentration of trypsinogen-2 in the sample exceeds the cut-off value for the test. A second blue line confirms correct performance of the test.

Performance of the test

Analytical Sensitivity

The analytical sensitivity, detection limit, of the Actim Pancreatitis test was evaluated using samples with different concentrations of human trypsinogen-2 on three different lots of Actim Pancreatitis test. The detection limit of the Actim Pancreatitis test is approximately 50 µg/l of human trypsinogen-2 in sample and the results remain positive at least to 100 000 µg/l.

Repeatability and reproducibility

Intra- (repeatability) and inter- (reproducibility) assay precision were evaluated using samples containing 0-100 000 µg/l of human trypsinogen-2 in urine. For repeatability the samples were tested on the same day using ten replicates at the concentrations of 0 – 100 µg/l and three replicates at the concentrations of 30 000 and 100 000 µg/l. The testing was done on three different lots of Actim Pancreatitis test, the total amount being 198 tests in one day. For reproducibility the samples were tested on the five consecutive workdays using two replicates per sample on three different lots of Actim Pancreatitis test, the total amount being 240 tests in five days. Repeatable and reproducible results were obtained.

Diagnostic performance

Actim Pancreatitis test has been evaluated in several clinical studies diagnosing acute pancreatitis with Actim Pancreatitis test. Selection of typical studies is shown in FIG 1 on the inner back cover.

Interference testing

The following interfering substances were tested with the Actim Pancreatitis test. No interference of the tested substances was observed with the performance of the Actim Pancreatitis test.

INTERFERING SUBSTANCE OR CONDITIONS	INTERFERENCE (NO/YES)
Glucose	No
Albumin	No
Hemoglobin	No
Acetylsalicylic acid	No
Ascorbic acid	No
Salicylic acid	No
Bilirubin	No
Biotin	No
pH 5.0-9.0	No

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Hedström J et al.

Urine trypsinogen-2 as a marker of acute pancreatitis. *Clin Chem* (1996) 42:685-690.

Jang T et al.

Point-of-care Urine Trypsinogen Testing for the Diagnosis of Pancreatitis. *Acad. Emerg. Med.* (2007) 14:29-34.

Kylänpää-Bäck M-L et al.

Reliable screening for acute pancreatitis with rapid urine trypsinogen-2 test strip. *Br J Surg* (2000), 87:49-52.

Kylänpää-Bäck M-L et al.

Comparison of urine trypsinogen-2 test strip with serum lipase in the diagnosis of acute pancreatitis. *Hepato-Gastroenterology* (2002) 49:1130-1134.

Kemppainen E et al.

Rapid measurement of urinary trypsinogen-2 as a screening test for acute pancreatitis. *N Engl J Med* (1997) 336: 1788-1793.

Actim[®]

Pancreatitis

Actim Oy
Klovinpellontie 3, FI-02180 Espoo, Finland
Tel. +358 9 547 680
www.medixbiochemica.com
medix@medixbiochemica.com

© Actim Oy. All rights reserved.
Actim is a trademark of Actim Oy.

AOACE32731-1 2019-04-24