
INSTRUCTIONS FOR USE





Actim[®] Calprotectin

EN	Instructions for use	3
CS	Návod k použití	6
DA	Brugsvejledning	9
DE	Gebrauchsanweisung	12
ES	Instrucciones de uso	15
FI	Käyttöohje	18
FR	Notice d'utilisation	21
IT	Istruzioni d'uso	24
NL	Gebruiksaanwijzing	27
NO	Bruksanvisning	30
PT	Instruções de utilização	33
SV	Bruksanvisning	36





Structure of dipstick

Specimen collection

Test procedure and results

Instructions for use

Numbers [1] - [12] refer to illustrations on inner cover.

Structure of dipstick

[1] Dip area [2] Result area [3] Test lines [4] Control line

Intended use

The Actim[®] Calprotectin test is a visually interpreted, semiquantitative immunochromatographic dipstick test, which is used for determination and monitoring of calprotectin concentrations in feces. The test is intended for professional use to help diagnose intestinal inflammation.

Kit components

The components of the Actim Calprotectin test are:

- Actim Calprotectin dipstick in a perforating device. The dipsticks (10) are packed in an aluminum container.
- Specimen Dilution Buffer (3 ml). This tris buffered solution contains bovine serum albumin (BSA), protease inhibitors and <0.1 % (w/v) sodium azide (NaN₃). One sampling stick is attached to the cap of each tube.

The kit 35032ETAC contains 2 aluminum tubes of dipsticks (35061ETAC) and 20 Specimen Dilution Buffer tubes (35071ETAC).

An instructions for use leaflet is included in all product formats. Toilet paper or a clean disposable container is needed for sampling (not provided with the test).

Storage

Store the dipsticks and the Specimen Dilution Buffer tubes at +2...+30 °C. Stored unopened, each component can be used until the expiry date marked on the component. After opening the aluminum container first time, the dipsticks can be stored for 4 months at +2...+30 °C. However, the expiry date must not be exceeded.

Specimen collection

The specimen is feces suspension that is prepared as soon as possible after sample collection. The original feces sample can be stored for 4 days at +2...+8 °C.

Unscrew the cap of the Specimen Dilution Buffer tube [5]. Collect the feces sample by using the sampling stick attached to the screw cap. There are two small slits in the head of each stick. Twist the stick randomly in several different

places in the feces sample so that both slits become filled with feces [6]. In case of liquid samples, dry the stick by wiping it with a piece of e.g. toilet paper before sticking it into the sample. The slits of the stick will be filled with the liquid feces.

Put the sampling stick into the Specimen Dilution Buffer tube by pushing the stick through the cone. Close the cap tightly [7] and shake the tube to suspend the feces in the buffer [8]. The specimen suspension is now ready for testing. The suspension can be stored for one week at +2...+25 °C (+2...+8 °C is recommended).

Test procedure and interpretation of the results

1. If stored refrigerated, allow the aluminum container of dipsticks and the buffer tube to reach room temperature before opening. Remove the required number of perforating devices (each containing a dipstick) from the container just prior to use and close the container immediately. Do not touch the yellow dip area at the lower ends of the dipsticks. Do not touch the part of the perforating device below the stopper, as this area will come into contact with the sample. Identifying marks may be written on the upper blue-grey parts of the dipsticks.
2. The specimen should be at room temperature upon testing. Use the sharp tip of the perforating device to puncture the perforation area in the bottom of the buffer tube **then push the device into the tube until the stopper prevents it going further [9]**.
3. Slowly invert the tube (not fully upside down) once and immediately place the tube on the table in an upright position [10]. **NOTE! It is important that the inverting movement is not too abrupt, a peaceful inversion is vital for the correct performance of the test.**
4. The result is available in 10 minutes [11]. **Do not pay attention to any lines appearing later than 10 minutes.**
5. **If, in addition to the blue control line in the upper part of the result area you see:**
 - **no other blue lines**, the calprotectin concentration is **<50 µg/g of feces**.
 - **one blue line in the lower part of the result area**, the calprotectin concentration is **50-200 µg/g of feces**.
 - **two blue lines**, the calprotectin concentration is **>200 µg/g of feces**.

If the control line does not appear, the test is invalid [12].

In case of a missing control line, please see the Notes section for troubleshooting.

The intensities of the blue lines will vary. When calprotectin concentration is low only one faint blue line in addition to the blue control line appears in 10 minutes. As the concentration of calprotectin increases, the intensity of the first line increases. A second blue line begins to appear when the calprotectin concentration exceeds 200 µg/g feces.

Disposal of the test device

The contact between the Specimen Dilution Buffer tube cap and the perforating device is so close that the test device can be disposed of directly after the result has been read. Do not remove the perforating device from the tube.

Limitations of the test

The test is intended for *in vitro* diagnostic use only.

Notes

- The test results are semiquantitative.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing.
- Do not touch the dip area of the dipstick or the area of the perforating device below the stopper.
- When inverting the sample tube to allow the dipstick to absorb sample, do not invert the tube too fast, or for too long. The test will not work properly if the amount absorbed is too small or too large.
- If the control line does not appear, the inversion of the sample tube may have been performed too abruptly. Take a new dipstick and repeat the inversion in a controlled way.
- Improper sampling may lead to false result.
- Calprotectin is not uniformly distributed in the feces sample. Therefore, it is important to dip the sampling stick randomly in several places in the feces sample.
- Appearance of a control line confirms correct performance of the test. If a control line does not appear the test is invalid, and should be repeated using another dipstick.

- If the test result cannot be interpreted clearly (e.g. if the lines are blotched or uneven) it is recommended that the test be repeated with another dipstick.
- At ten minutes the appearance of any faint-to-dark blue line along with a control line indicates a positive result. However, do not pay attention to any lines appearing after 10 minutes.
- If only the control line is visible, the test result should be interpreted as negative only after 10 minutes have elapsed.
- If the color of the test line is not blue, the result should be interpreted as invalid and the test should be repeated with a new feces sample.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

Principle of the test

Calprotectin is a degradation product from neutrophil granulocytes from mucosal layer of the colon. It is released in excess into the intestines when colonic inflammation is present. Fecal calprotectin testing can be used as an aid to distinguish between inflammatory bowel diseases, such as Crohn's disease and ulcerative colitis, and non-inflammatory bowel diseases, such as irritable bowel syndrome. It can also be used to determine the disease activity and monitor the response to treatment in patients with ulcerative colitis or Crohn's disease.

The test is based on immunochromatography. It involves monoclonal antibodies to human calprotectin. Antibodies are bound to blue latex particles (the detecting label). In the carrier membrane there are two calprotectin specific antibody zones, to which the blue latex particles will bind if the sample contains calprotectin.

When the feces suspension comes into contact with the dip area of the dipstick, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains calprotectin, it binds to the antibodies attached to the latex particles. The particles are carried by the liquid flow and, if calprotectin are bound to them, they bind to the corresponding catching antibody zones. If sample contains >200 µg calprotectin/g feces two blue lines will appear. A blue control line confirms correct performance of the test.

Performance of the test

Analytical Sensitivity

The analytical sensitivity, detection limit, of the Actim Calprotectin test was evaluated by two operators using samples with different concentrations of human calprotectin corresponding to the 12.5 – 10 000 µg of calprotectin/g of feces on three different lots of Actim Calprotectin test.

The detection limit of Actim Calprotectin is 50 µg/g of feces. At this point 97 % of the results were positive. The second blue line begins to appear when feces suspension contains calprotectin > 200 µg/g of feces. At this point 90 % of the results were positive. Hook effect has been tested with sample containing up to 10 000 µg of human calprotectin/g of feces giving consistently the positive result.

Repeatability and reproducibility

Intra- (repeatability) and inter- (reproducibility) assay precision were evaluated using samples containing 0 – 200 µg of human calprotectin/g of feces. For repeatability the samples were tested by two operators on the same day using 10 replicates per sample on three different lots of Actim Calprotectin, the total amount being 120 tests/Lot (360/3 lots) in one day. The agreement for the samples containing 50 µg of human calprotectin/g of feces was 97 % and for the sample, containing 200 µg of human calprotectin/g of feces was 97 %. For reproducibility the samples were tested on the five consecutive workdays using three replicates per sample on three different lots of Actim Calprotectin, the total amount being 180 tests/Lot (540/3 lots) in five days. The agreement for the samples containing 50 µg of human calprotectin/g of feces was 98 % and for the sample, containing 200 µg of human calprotectin/g of feces was 98 %.

Analytical specificity, Cross-reactivity

The Actim Calprotectin test was tested with cross-reactivity's of human and bovine transferrin and hemoglobin. No cross-reactivity was observed. The Actim Calprotectin test is specific for the human calprotectin.

Method comparison

The function of Actim Calprotectin test was evaluated by testing 100 patient samples with Actim Calprotectin and measuring the calprotectin concentration from the same fecal suspensions by Bühlmann fCAL ELISA. The evaluation of the results showed that the Actim Calprotectin results are comparable to the quantitative measurement. The results are shown in FIG 1 on the inner back cover.

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© Actim Oy. All rights reserved.
Actim is a trademark of Actim Oy.
Actim Oy