External controls for Actim Combi® Hemoglobin Transferrin test, Instructions for use

INTENDED USE

Actim Combi® Hemoglobin Transferrin Controls are intended to be used with the Actim Combi® Hemoglobin Transferrin test (30322ETAC) for external quality control. The controls may also be used to demonstrate negative results and weak and strong positive results.

KIT COMPONENTS

Actim Combi Hemoglobin Transferrin Controls kit (30310ETAC) consists of one of each of the following components:

- Actim Combi Hemoglobin Transferrin Negative Control (30390ETAC)
- Actim Combi Hemoglobin Transferrin Positive Control, Low (30391ETAC)
- Actim Combi Hemoglobin Transferrin Positive Control, High (30392ETAC)
- Actim Reconstitution Solution (30889ETAC)
- Instructions for use

A pipette capable to measure 500 µl of solution is needed for accurate dispensing of Reconstitution Solution (not provided with the kit). For dilution of the controls, a pipette capable to measure 100 µl is also needed.

STORAGE

Store the components at +2...+8 °C. Stored unopened, each component can be used until the expiry date marked on the component. After reconstitution, the controls should be used within 7 days. Alternatively, the reconstituted controls can also be stored frozen at -20 °C for 12 months. Do not freeze after thawing.

CONTROL MATERIAL

The Positive Control consist of purified human hemoglobin and transferrin in buffered solution containing protein as stabilizer and <0.1 % sodium azide (NaN₃). The Negative Control consists of the same matrix without added antigen.

LIMITATIONS OF THE TEST

- The controls are to be used with the Actim Combi Hemoglobin Transferrin test (manufactured by Actim Oy) only.
- The controls are intended for professional *in vitro* use only.

TEST PROCEDURE AND INTERPRETATION OF THE RESULTS

Reconstitution:

1. Unscrew the cap and open the vials by carefully removing the grey stopper.
2. Reconstitute the controls by adding 500 µl of Reconstitution Solution to the vials. When pipetting, avoid the risk of carryover by changing the pipette or its tip after each addition.
3. Close the stopper and allow the vial to reconstitute for at least 15 minutes.
4. Mix the solutions.

Testing:

1. Dilute each control by adding 100 µl of the control to the Specimen Dilution Buffer (30371ETAC). When pipetting, avoid the risk of carryover by changing the pipette or its tip after each addition.
2. Testing is performed by carefully following the instructions for use of the Actim Combi Hemoglobin Transferrin test.

The test results are interpreted as instructed in the paragraph "Test procedure and interpretation of results" of the Actim Combi Hemoglobin Transferrin test:

- The negative control gives a negative test result.
- The low positive control gives a weak positive result with hemoglobin line and transferrin line.
- The high positive control gives a strong positive result with hemoglobin line and transferrin line.

NOTES

- Before use, ensure that the package and vials are intact. Do not use a control if the vial is damaged.
- Do not use the controls before they have been completely reconstituted.
- When testing, the controls should be at room temperature to ensure optimal test results with Actim Combi Hemoglobin Transferrin test.
- Carefully follow the instructions for use of the Actim Combi Hemoglobin Transferrin test and Actim Combi Hemoglobin Transferrin Controls to ensure correct results.
- If the results of the controls are not appropriate, repeat the procedure with another dipstick.
- The frequency of quality control testing is to be determined by each laboratory independently, and will depend on factors such as patient sample frequency and/or local regulations.
- Plasma from each donor has been tested and found negative for antibody to HIV-1/ HIV-2, antibody to HCV, for HBSAG, and syphilis. Nevertheless, such tests are unable to prove the complete absence of viruses, and the controls should therefore be handled taking appropriate precautions.
- Like all biological materials, the controls should be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

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