
INSTRUCTIONS FOR USE




Actim Combi[®] Hemoglobin. Transferrin.

EN	Instructions for use	3
CS	Návod k použití	6
DA	Brugsvejledning	9
DE	Gebrauchsanweisung	12
ES	Instrucciones de uso	15
FI	Käyttöohje	18
FR	Notice d'utilisation	21
IT	Istruzioni d'uso	24
NL	Gebruiksaanwijzing	27
NO	Bruksanvisning	30
PT	Instruções de utilização	33
SV	Bruksanvisning	36





Structure of dipstick

Specimen collection

Test procedure and results

Instructions for use

Numbers [1] - [12] refer to illustrations on inner cover.

Structure of dipstick

[1] Dip area [2] Result area [3] Test lines [4] Control line

Intended use

The Actim Combi[®] Hemoglobin Transferrin test is a visually interpreted, qualitative immunochromatographic dipstick test for detecting occult blood in feces. The Actim Combi Hemoglobin Transferrin test detects human hemoglobin and human transferrin in a feces suspension. The test is intended for professional use to help diagnose gastrointestinal bleeding.

Kit components

The components of the Actim Combi Hemoglobin Transferrin test are:

- Actim Combi Hemoglobin Transferrin dipstick in a perforating device. The dipsticks (10) are packed in an aluminum container.
- Specimen Dilution Buffer (3 ml). This phosphatebuffered solution contains bovine serum albumin (BSA), protease inhibitors and <0.1 % (w/v) sodium azide (NaN₃). One sampling stick is attached to the cap of each tube.

The kit 30322ETAC contains 2 aluminum tubes of dipsticks (30351ETAC) and 20 Specimen Dilution Buffer tubes (30371ETAC).

An instructions for use leaflet is included in all product formats. Toilet paper or a clean disposable container is needed for sampling (not provided with the test).

Storage

Store the dipsticks and the Specimen Dilution Buffer tubes at +2...+30 °C. Stored unopened, each component can be used until the expiry date marked on the component. After opening the aluminum container first time, the dipsticks can be stored for 4 months at +2...+30 °C. However, the expiry date must not be exceeded.

Specimen collection

The specimen is feces suspension that is prepared as soon as possible after sample collection. The original feces sample can be stored for 3 days at +2...+8 °C.

Unscrew the cap of the Specimen Dilution Buffer tube [5]. Collect the feces sample by using the sampling stick attached to the screw cap. There are two small slits in the head of each stick. Twist the stick randomly in several different places in the feces sample so that both slits become filled with feces [6]. In case of liquid samples, dry the stick by wiping it with a piece of e.g. toilet paper before sticking it into the sample. The slits of the stick will be filled with the liquid feces.

Put the sampling stick into the Specimen Dilution Buffer tube by pushing the stick through the cone. Close the cap tightly [7] and shake the tube to suspend the feces in the buffer [8]. The specimen suspension is now ready for testing. The suspension can be stored for one week at +2...+25 °C (+2...+8 °C is recommended).

Test procedure and interpretation of the results

1. If stored refrigerated, allow the aluminum container of dipsticks and the buffer tube to reach room temperature before opening. Remove the required number of perforating devices (each containing a dipstick) from the container just prior to use and close the container immediately. Do not touch the yellow dip area at the lower ends of the dipsticks. Do not touch the part of the perforating device below the stopper, as this area will come into contact with the sample. Identifying marks may be written on the upper red parts of the dipsticks.
2. The specimen should be at room temperature upon testing. Use the sharp tip of the perforating device to puncture the perforation area in the bottom of the buffer tube **then push the device into the tube until the stopper prevents it going further [9].**
3. Slowly invert the tube (not fully upside down) once and immediately place the tube on the table in an upright position [10]. **NOTE! It is important that the inverting movement is not too abrupt, a peaceful inversion is vital for the correct performance of the test.**
4. The result can be interpreted as positive as soon as one blue line or one red and or both blue and red lines, and one black line become visible in the result area. The negative result should be read at **10 minutes** after the liquid front has reached the lower part of the result area [11]. **Do not pay attention to any lines appearing later than 10 minutes.**

5. If, in addition to the black control line you see:

- **one blue line** in the **lower part** of the result area, the result is **hemoglobin positive**.
- **one red line** in the **middle part** of the result area, the result is **transferrin positive**.
- **one blue line and one red line** in the result area, the result is **hemoglobin and transferrin positive**.
- **no blue or red lines** in the result area, the result is **negative**.

If the control line does not appear, the test is **invalid [12]**.

In case of a missing control line, please see the Notes section for troubleshooting.

Disposal of the test device

The contact between the Specimen Dilution Buffer tube cap and the perforating device is so close that the test device can be disposed of directly after the result has been read. Do not remove the perforating device from the tube.

Limitations of the test

The test is intended for *in vitro* diagnostic use only.

Notes

- The test results are qualitative. No quantitative interpretation should be made based on the test results.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing.
- Take care that any part of the test device that will be in contact with the sample tube do not become contaminated by external blood, e.g. from a cut finger. Do not touch the dip area of the dipstick or the area of the perforating device below the stopper.
- When inverting the sample tube to allow the dipstick to absorb sample, do not invert the tube too fast, or for too long. The test will not work properly if the amount absorbed is too small or too large.
- If the control line does not appear, the inversion of the sample tube may have been performed too abruptly. Take a new dipstick and repeat the inversion in a controlled way.
- Improper sampling may lead to false result.
- Occult blood is not uniformly distributed in the feces sample. Therefore, it is important to dip the sampling stick randomly in several places in the feces sample.
- Appearance of a control line confirms correct performance of the test. If a control line does not appear the test is invalid, and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly (e.g. if the lines are blotched or uneven) it is recommended that the test be repeated with another dipstick.
- At ten minutes the appearance of any faint-to-dark blue and/or red test line along with a control line indicates a positive result. However, do not pay attention to any lines appearing after 10 minutes.
- If only the control line is visible, the test result should be interpreted as negative only after 10 minutes have elapsed.
- If the color of the test line is not blue (hemoglobin) or red (transferrin), the result should be interpreted as invalid and the test should be repeated with a new feces sample.
- Patients with the following conditions should not be considered for testing, as these conditions may interfere with test results: bleeding hemorrhoids, constipation bleeding, menstrual bleeding and urinary bleeding. However, these patients may be considered for testing after such bleeding ceases.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

Principle of the test

Feces do not normally contain measurable amounts of blood. The presence of hemoglobin and transferrin in feces implies that there may be bleeding in the gastrointestinal tract.

The test is based on immunochromatography. It involves antibodies to human hemoglobin and to human transferrin. Antibodies are bound to blue (hemoglobin) or red (transferrin) latex particles (the detecting label). Antibodies against both types are immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result.

When the feces suspension comes into contact with the dip area of the dipstick, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains hemoglobin and /or transferrin, they bind to the antibodies attached to the latex particles. The particles are carried by the liquid flow and, if hemoglobin and / or transferrin are bound to them, they bind to the corresponding catching antibodies. A blue line will appear in the result area if the concentration of hemoglobin in the sample exceeds the detection limit of the test. A red line will appear in the result area if the concentration of transferrin in the sample exceeds the detection limit of the test. A black line (control line) confirms correct performance of the test.

Performance of the test

Analytical Sensitivity

The analytical sensitivity, detection limit, of the Actim Combi Hemoglobin Transferrin test was evaluated by two operators using samples with different concentrations ranging of human hemoglobin from 10 – 250 000 µg/l and human transferrin from 1 – 10 000 µg/l on three different lots of Actim Combi Hemoglobin Transferrin test.

The detection limit of the Actim Combi Hemoglobin Transferrin test is approximately 50 µg/l of human hemoglobin in the sample and the results remain positive at least to 250 000 µg/l. A feces suspension with 50 µg of hemoglobin/l corresponds to 25 µg of hemoglobin/g of feces.

The detection limit of the Actim Combi Hemoglobin Transferrin test is approximately 5 µg/l of human transferrin in the sample and the results remain positive at least to 10 000 µg/l. A feces suspension with 5 µg/l of transferrin corresponds to 1.5 µg of transferrin /g of feces. Transferrin 2.5 µg/l can be positive or negative which corresponds to 0.75 µg/g of feces.

Repeatability and reproducibility

Intra- (repeatability) and inter- (reproducibility) assay precision were evaluated using samples containing 0 – 200 µg/l of human hemoglobin and transferrin in samples. For repeatability the samples were tested by two operators on the same day using 10 replicates per sample on three different lots of Actim Combi Hemoglobin Transferrin, the total amount being 140 tests/Lot (420/3 lots) in one day. For reproducibility the samples were tested on the five consecutive workdays using three replicates per sample on three different lots of Actim Combi Hemoglobin Transferrin, the total amount being 210 tests/Lot (630/3 lots) in five days. Repeatable and reproducible results were obtained.

Analytical specificity, Cross-reactivity

The Actim Combi Hemoglobin Transferrin test was tested with following cross-reactivity's: Cross-reactivity of human transferrin in hemoglobin-line, cross-reactivity of human hemoglobin in transferrin-line, cross-reactivity of Bovine hemoglobin, Porcine hemoglobin, Bovine transferrin, Porcine transferrin and Human Lactoferrin. No cross-reactivity was observed. The Actim Combi Hemoglobin Transferrin test is specific for the human hemoglobin and transferrin.

Method comparison

The function of Actim Combi Hemoglobin Transferrin test was compared to the immunochromatographic test of CerTest FOB-Transferrin. The agreement with the two methods was 97.1 % (hemoglobin) and 93.3 % (transferrin). The results are shown in FIG 1 on the inner back cover.

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FIG 1.


COMPARISON RESULTS OF THE HEMOGLOBIN LINE:

Actim Combi® Hemoglobin Transferrin	Reference: CerTest FOB-Transferrin		
	POS	NEG	TOTAL
POS	26	3	29
NEG	0	76	76
OVERALL AGREEMENT 97.1%			n=105

COMPARISON RESULTS OF THE TRANSFERRIN LINE

Actim Combi® Hemoglobin Transferrin	Reference: CerTest FOB-Transferrin		
	POS	NEG	TOTAL
POS	35	3	38
NEG	4	63	67
OVERALL AGREEMENT 93.3%			n=105

Explanation of symbols • Vysvětlení symbolů • Symbolforklaring • Interpretation der Symbole • Significado de los símbolos • Symbolien selitykset • Explication des symboles • Spiegazione dei simboli • Interpretatie van de symbolen • Symbolforklaring • Explicação dos símbolos • Teckenförklaringar

YYYY-MM-DD

Use By
Použitelné do
Anvendes inden
Verwendbar bis
Fecha de caducidad
Viimeinen käyttöpäivä
Utiliser jusque
Utilizzare entro
Houdbaar tot
Brukes før
Válido até
Använd före

Temperature limitation
Teplotní rozmezí od do
Temperaturbegrensning
Temperaturbereich
Limite de temperatura
Lämpötilarajoitus
Limites de température
Limiti di temperatura
Temperatuurlimiet
Temperaturbegrensning
Limite da temperatura
Temperaturbegrensning

Manufacturer
Výrobce
Producent
Hersteller
Fabricante
Valmistaja
Fabricant
Fabbriicante
Fabrikant
Produsent
Fabricante
Tillverkare

LOT

Batch code
Číslo šarže
Batchnummer
Chargenbezeichnung
Código de lote
Erä
Code du lot
Codice del lotto
Lot nummer
Lotnummer
Número de lote
Lot nummer

REF

Catalogue number
Katalogové číslo
Katalognummer
Bestellnummer
Número de catálogo
Tuotenumero
Référéncia du catalogue
Numero di catalogo
Catalogusnummer
Katalognummer
Código do produto
Katalognummer

Single use
Jednorázové použití
Engangsbrug
Nicht wiederverwendbar
No reutilizar
Kertakäyttöinen
Ne pas réutiliser
Non riutilizzare
Niet opnieuw gebruiken
Ikke til gjenbruk
De uso único
Återanvänds ej

Contains sufficient for <n> tests
Lze použít pro <n> testů
Indeholder tilstrækkeligt til "n" test
Ausreichend für "n" Tests
Contenido suficiente para "n" ensayos
Sisältö on riittävä "n" testiä varten
Contenu suffisant pour "n" tests
Contenuto sufficiente per "n" saggi
Inhoud voldoende voor "n" testen
Inneholder nok til "n" tester
Conteúdo suficiente para <n> exames
Räcker till "n" antal tester

DILSPE

Specimen dilution buffer
Pufř k ředění vzorku
Fortyndingsbuffer
Probenverdünnungspuffer
Tampón de dilución de Muestra
Näytteenlaimennuspuskuri
Tampon de dilution
Tampone di diluzione
Buffer om de monsters te verdunnen
Fortynningsbuffer
Tampão para diluição da amostra
Provutspädningslösning

STICK

Dipstick
Testovací proužek
Teststrimmel
Teststreifen
Tira reactiva
Testitikk
Bandelette
Dipstick
Dipstick
Teststrimmel
Tira reactiva
Teststicka

Consult instructions for use
Viz návod k použití
Se brugsanvisningen
Gebrauchsanweisung zurate ziehen
Consultar las instrucciones de uso
Katso käyttöohjeet
Consulter la notice d'utilisation
Consultare le istruzioni d'uso
Raadpleeg de gebruiksinstructies
Se bruksanvisningen
Consultar as instruções de utilização
Läs i bruksanvisningen

IVD

In Vitro Diagnostic Medical Device
In Vitro diagnostický zdravotnický prostředek
Medicinsk udstyr til *in vitro*-diagnostik
In vitro Diagnostikum
Producto sanitario para diagnóstico *in vitro*
In vitro -diagnostiikkaan tarkoitettu
lääkinnällinen laite
Dispositif médical de diagnostic *in vitro*
Dispositivo medico-diagnostico *in vitro*
Medisch hulpmiddel voor *in-vitro*diagnostiek
Medisinsk utstyr for *in vitro* diagnostikk
Dispositivos médicos de diagnóstico *in vitro*
Medicintekniska produkter för *in vitro* diagnostik

CONT NaN₃ <0.1%

Contains <0.1 % NaN₃
Obsahuje <0,1 % NaN₃
Indeholder <0,1 % NaN₃
Enthält <0,1 % NaN₃
Contiene <0,1 % NaN₃
Sisältää <0,1 % NaN₃
Contient <0,1 % NaN₃
Contiene <0,1 % NaN₃
Bevat <0,1 % NaN₃
Inneholder <0,1 % NaN₃
Contém <0,1 % NaN₃
Innehåller <0,1 % NaN₃

Actim Combi[®]

Hemoglobin. Transferrin.

Actim Oy
Klovinpellontie 3, FI-02180 Espoo, Finland
Tel. +358 9 547 680
www.actimtest.com
actim@actimtest.com

© Actim Oy. All rights reserved.
Actim is a trademark of Actim Oy.

AOACE30322-2 2020-06-30