## INSTRUCTIONS FOR USE


# Actim<sup>®</sup> Influenza A&B

<b>EN</b>	Instructions for use	<b>3</b>
<b>DE</b>	Gebrauchsanweisung	<b>6</b>
<b>ES</b>	Instrucciones de uso	<b>10</b>
<b>FI</b>	Käyttöohje	<b>13</b>
<b>FR</b>	Notice d'utilisation	<b>16</b>
<b>IT</b>	Istruzioni d'uso	<b>20</b>
<b>PL</b>	Instrukcja użytkowania	<b>23</b>
<b>PT</b>	Instruções de utilização	<b>26</b>
<b>SV</b>	Bruksanvisning	<b>29</b>

## Structure of dipstick

## Test procedure and results

# Instructions for use

Numbers ① - ⑨ refer to illustrations on inner cover.

## Structure of dipstick

① Dip area ② Result area ③ Test lines ④ Control line

## Intended use

The Actim® Influenza A&B test is a visually interpreted, qualitative immunochromatographic rapid test for detection of influenza type A and type B virus antigens in respiratory samples. The test is intended for professional use to help diagnose acute influenza infections.

## Kit components

The components of the Actim Influenza A&B test kits 32836ETAC and 32837ETAC are:

- 20 individually packed dipsticks in sealed aluminum foil pouches with desiccant.
- One dropper bottle of Swab Extraction Buffer (10 ml). This solution contains detergents, bovine serum albumin (BSA) and preservatives.
- One dropper bottle of Extraction Buffer (2.3 ml). This solution contains detergents, bovine serum albumin (BSA) and preservatives.
- 20 swabs for swab sample collection.
- 20 test tubes for sample extraction.
- Instructions for use and visual test procedure.

In addition Actim Influenza A&B test kit 32836ETAC contains external influenza control swabs: one control swab each of Positive Control A and Positive Control B.

Pipettes may be needed for sampling (not provided in the test kit).

## Storage

Store the test kit at +2...+30 °C. Each component can be used until the expiry date marked on the component. Use the dipsticks and controls shortly after their removal from the aluminum foil pouch, because moisture damages them.

## Specimen collection and handling

The specimen is a respiratory sample. Collect a nasal swab sample, a nasopharyngeal swab sample, a nasal aspirate or a nasal wash sample. The samples should be tested as soon as possible after collection.

If transport of the samples is required, the swab sample can be placed in a test tube and stored dry for 6 hours at room temperature before testing. Alternatively, the following transport media have been tested and shown not to affect the function of the test: UTM-RT Transport medium, Hank's Balanced Salt Solution, Saline Solution, Virocult and Dulbecco's Modified Eagle Medium with 1 % FCS supplemented with penicillin and gentamycin. The transported samples should be tested within two days.

## Nasal/nasopharyngeal swab procedure:

1. Slowly invert the dropper bottle of swab extraction buffer upside down and dispense 12 drops into a test tube (= extraction tube) ⑤.
2. Place the swab into the test tube. The specimen is extracted from the swab by swirling the latter vigorously for approximately 10 seconds ⑥. Press the swab against the wall of the tube and remove it from the tube.

## Nasal wash/aspirate or sample in transport media procedure:

3. Dispense 300 µl of nasal wash or nasal aspirate solution into a test tube (= extraction tube).
4. Slowly invert the dropper bottle of extraction buffer upside down and dispense 3 drops into the same test tube ⑤. Mix the solutions well ⑥.

Alternatively, if the aspirate volume is very small, a swab can be used to transfer the specimen to a test tube containing 12 drops of swab extraction buffer. Proceed according to the nasal/nasopharyngeal swab procedure. Continue to the test procedure.

## Test procedure and interpretation of the results

1. Allow the aluminum foil pouch to reach room temperature if stored refrigerated. Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper blue part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch.
2. Place the yellow dip area [the lower end of the dipstick] into the test tube with extracted sample and leave it there until reading ⑦.

3. Read the result **at 10 minutes**. A result can be interpreted as positive as soon as one or two blue and one black line become visible in the result area. If only one black line appears within 10 minutes, the result is interpreted as negative 8.

Appearance of a black control line confirms that the test has been performed properly.

If, in addition to the black control line you see:

- **one blue line** in the **lower part** of the result area, the result is **influenza A positive**.
- **one blue line** in the **middle part** of the result area, the result is **influenza B positive**.
- **two blue lines** in the result area, the result is **influenza A and B positive**.
- **no blue lines** in the result area, the result is **influenza A and B negative**.

If the control line does not appear the test is invalid 9.

#### External quality control (only with the kit 32836ETAC)

Perform according to the Nasal/nasopharyngeal swab procedure.

1. Slowly invert the dropper bottle of swab extraction buffer upside down and dispense 12 drops into a test tube (= extraction tube).
2. Place the control swab into the test tube. The control is extracted from the control swab by swirling the latter vigorously for approximately 10 seconds. Press the swab against the wall of the tube and remove it from the tube.
3. Proceed to the test procedure.

#### Limitations of the test

The test is intended for *in vitro* diagnostic use only.

#### Notes

- No quantitative interpretation should be made based on the test results.
- Use only the swab provided with the kit.
- Nasal aspirate samples may be very mucous. Such samples should be pipetted carefully to ensure that the pipette absorbs 300 µl of sample. Very mucous samples are occasionally not absorbed by the dipstick. Such specimens may be diluted 1:2 with swab extraction buffer, and retested.
- The amount of virus in the extracted sample may remain lower than usual if

the sample needs to be diluted or transferred with the swab or extracted into a viral transport media. Nasal wash samples may be dilute and contain a low amount of virus. In these cases the results may be less strongly positive than usual, or even remain negative.

- A negative test result may occur if the sample was collected or stored improperly.
- Positive test results do not rule out co-infections with other pathogens.
- The test detects both viable and non-viable influenza virus.
- Do not use the dipstick, control swab or swab if their foil pouches or seals of the pouches are not intact.
- Care must be taken when placing the dipstick in the sample. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the test.
- Do not use a dipstick if you notice a blue or black coloring in the result area before testing.
- When using the dropper bottles, ensure that the bottle is completely inverted upside down before dropping. This optimizes the drop size.
- To avoid foaming, do not shake the extraction buffer and swab extraction buffer dropper bottles.
- If a control line does not appear the test is invalid, and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly it is recommended that the test be repeated with another dipstick.
- The appearance of any faint-to-dark blue test line along with a black control line indicates a positive result.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

#### Principle of the test

The test includes the extraction and detection of influenza type A and B virus antigens. The test is based on immunochromatography. It utilizes monoclonal

antibodies to influenza type A and type B virus antigens. Antibodies against both types are bound to blue latex particles (the detecting label). Antibodies against both types are immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result.

When the dip area of the dipstick is placed in an extracted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains influenza type A or B, the virus antigen binds to either of the antibodies labeled with latex particles. The particles are carried by the liquid flow and, if influenza virus antigen is bound to them, they bind to the corresponding catching antibody on the membrane. A blue line (test line) will appear in the lower part of the result area if the concentration of influenza type A in the sample exceeds the detection limit of the test. A blue line (test line) will appear in the middle of the result area if the concentration of influenza type B in the sample exceeds the detection limit of the test. Appearance of the black control line confirms correct performance of the test.

## Performance of the test

### Analytical Sensitivity

The analytical sensitivity of Actim Influenza A&B test was assessed using different influenza viral strains. The detection limit of influenza virus is  $1.0 \times 10^4$  TCID<sub>50</sub>/ml for influenza A/New Caledonia/20/99 (H1N1),  $1.0 \times 10^4$  TCID<sub>50</sub>/ml for influenza A/Panama/2007/99 (H3N2) and  $2.5 \times 10^4$  TCID<sub>50</sub>/ml for influenza B/Yamanashi/166/98. The test has further been shown to detect cultured avian influenza A subtypes (H5N1, H5N3, H7N3, H7N9, H9N2).

### Analytical Specificity

Actim Influenza A&B test was evaluated with the microorganisms listed below. A negative test result was obtained with all microorganisms. Test results indicate that Actim Influenza A&B test is specific for influenza A and influenza B.

### Viruses/Bacteria ( $\geq 10,000$ TCID<sub>50</sub>/ml)

- Echovirus types 2, 3, 6, 9, 11, 25, 30
- Adenovirus types 1, 2, 3, 4, 5, 6, 7, 10, 18
- Parainfluenza types 1, 2, 3
- Coxsackievirus types A9, B5, B6
- Rhinovirus types 1A, 2, 13, 15, 37
- Respiratory syncytial virus (RSV) types A and B
- Cytomegalovirus AD169
- Herpes simplex virus types 1 and 2
- *Mycoplasma pneumoniae* ( $>10,000$  CFU)

### Bacteria/Yeast ( $\geq 10^5$ cfu/ml)

- *Bordetella pertussis*
- *Bacteroides fragilis*
- *Neisseria gonorrhoeae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Candida albicans*
- *Streptococcus viridians*
- *Streptococcus pneumoniae*
- *Listeria monocytogenes*
- *Streptococcus pyogenes*
- *Streptococcus agalactiae*
- *Streptococcus ssp. Group C*
- *Streptococcus ssp. Group G*
- *Moraxella catarrhalis*
- *Enterococcus faecalis*
- *Staphylococcus Aureus ATCC 25923*
- *Staphylococcus epidermidis*
- *Klebsiella pneumoniae*
- *Pseudomonas aeruginosa*
- *Escherichia coli*
- *Proteus mirabilis*
- *Staphylococcus aureus ATCC 29213*
- *Staphylococcus aureus ATCC 43300*

### Repeatability and Reproducibility

A panel of specimens consisting of samples of influenza A (low positive and high positive), samples of influenza B (low positive and high positive) and a negative sample was evaluated for intra- (repeatability) and inter- (reproducibility) assay precision. For repeatability the samples were tested with 10 replicates during the same day using one lot of the Actim Influenza A&B test. For reproducibility, the samples were tested on 7 days using one lot of the Actim Influenza A&B test. One hundred percent (100%) agreement with expected test results was obtained.

### Interfering substances

Nasal drops (30 %) and ibuprofen (15 mg/ml) were tested with Actim Influenza A&B test. No interference of these substances was observed with the test performance. Whole blood was tested with Actim Influenza A&B test. No interference of whole blood (2%) was observed with the test performance.

### Diagnostic performance

Actim Influenza A&B test was evaluated at Institute of Medical & Veterinary Science, Adelaide, Australia and at University of Turku, Department of Virology, Finland. The results are shown in FIG 1. and FIG 2. on the inner back cover.

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\*Due to limited number of positive influenza B infections in the study population, it is impossible to define the clinical performance.

FIG 2.

<b>INFLUENZA A</b>				<b>INFLUENZA B</b>			
<b>Actim Influenza A&amp;B</b>	<b>REFERENCE METHOD TR-FIA</b>			<b>Actim Influenza A&amp;B</b>	<b>REFERENCE METHOD TR-FIA</b>		
	POS	NEG	TOTAL		POS	NEG	TOTAL
POS	18	0	18	POS	19	0	19
NEG	3	68	71	NEG	2	68	70
TOTAL	21	68	89	TOTAL	21	68	89
		<b>95 % CI</b>				<b>95 % CI</b>	
SENS	86 %	63–96 %		SENS	90 %	68–98 %	
SPEC	100 %	93–100 %		SPEC	100 %	93–100 %	
PPV	100 %	78–100 %		PPV	100 %	79–100 %	
NPV	96 %	87–99 %		NPV	97 %	89–100 %	

**Explanation of symbols** • Interpretation der Symbole • Significado de los símbolos • Symbolien selitykset • Explication des symboles • Spiegazione dei simboli •  
Wyjaśnienie symboli • Explicação dos símbolos • Teckenförklaringar

YYYY-MM-DD

Use By  
Verwendbar bis  
Fecha de caducidad  
Viimeinen käyttöpäivä  
Utiliser jusque  
Utilizzare entro  
Użyć przed  
Válido até  
Använd före

Temperature limitation  
Temperaturbereich  
Limite de temperatura  
Lämpötilarajoitus  
Limites de température  
Limiti di temperatura  
Przestrzegać zakresu temperatury  
Limite da temperatura  
Temperaturbegränsning

Manufacturer  
Hersteller  
Fabricante  
Valmistaja  
Fabricant  
Fabricante  
Produttore  
Fabricante  
Tillverkare

**LOT**

Batch code  
Chargenbezeichnung  
Codigo de lote  
Erä  
Code du lot  
Codice del lotto  
Kod partii  
Número de lote  
Lot number

**REF**

Catalogue number  
Bestellnummer  
Número de catálogo  
Tuotenumero  
Référence du catalogue  
Numero di catalogo  
Numer katalogowy  
Código do produto  
Katalognummer

Swab  
Tupfer  
Hisopo  
Näytteenottotikku  
Écouvillon  
Tampone  
Wymaz  
Zaragatoa  
Provtagningsspinne

**SWAB EXTR**

Extraction buffer  
Extraktionspuffer  
Tampón de extracción  
Uuttopuskuri  
Tampon d'extraction  
Tampone di estrazione  
Bufor ekstrakcyjny  
Tampão de Extração  
Extraktionsbuffert

**CONTROL +**

Positive Control  
Positive Kontrolle  
Control positivo  
Positiivinen kontrolli  
Contrôle positif  
Controllo positivo  
Kontrola pozytywna  
Controlo de Influenza A  
Positiv kontroll

**CONTROL A**

Influenza A Control  
Influenza A Kontrolle  
Control de Influenza A  
Influenssa A kontrolli  
Contrôle pour la grippe de type A  
Controllo dell'influenza A  
Kontrola pozytywna grypy typu A  
Controlo de Influenza A  
Influenza A kontroll

**CONTROL B**

Influenza B Control  
Influenza B Kontrolle  
Control de Influenza B  
Influenssa B kontrolli  
Contrôle pour la grippe de type B  
Controllo dell'influenza B  
Kontrola pozytywna grypy typu B  
Controlo de Influenza B  
Influenza B kontroll

**BUF EXTR**

Extraction buffer  
Extraktionspuffer  
Tampón de extracción  
Uuttopuskuri  
Tampon d'extraction  
Tampone di estrazione  
Bufor ekstrakcyjny  
Tampão de Extração  
Extraktionsbuffert

Consult instructions for use  
Gebrauchsanweisung zurate ziehen  
Consultar las instrucciones de uso  
Katso käyttöohjeet  
Consulter la notice d'utilisation  
Consultare le istruzioni d'uso  
Należy zapoznać się z instrukcją użytkowania  
Consultar as instruções de utilização  
Läs i bruksanvisningen

Contains sufficient for <n> tests  
Ausreichend für "n" Tests  
Contenido suficiente para "n" ensayos  
Sisältö on riittävä "n" testiä varten  
Contenu suffisant pour "n" tests  
Contenuto sufficiente per "n" saggi  
Wystarczy na wykonanie <n> testów  
Conteúdo suficiente para <n> exames  
Räcker till "n" antal tester

**STERILE EO**

Sterilized using ethylene oxide  
Sterilisation mit Ethylenoxid  
Esterilizado con óxido de etileno  
Steriloitu etyleenioksidilla  
Stérilisé à l'oxyde d'éthylène  
Sterilizzato con ossido di etilene  
Sterylizowane za pomocą tlenku etylenu  
Esterilizado com óxido de etileno  
Steriliserad med etylenoxid

**IVD**

*In Vitro* Diagnostic Medical Device  
*In vitro* Diagnostikum  
Producto sanitario para diagnóstico *in vitro*  
*In vitro* -diagnostiikkaan tarkoitettu lääkinällinen laite  
Dispositif médical de diagnostic *in vitro*  
Dispositivo medico-diagnostico *in vitro*  
Wyrób do diagnostyki *in vitro*  
Dispositivos médicos de diagnóstico *in vitro*  
Medicintekniska produkter för *in vitro* diagnostik

Single use  
Nicht wiederverwendbar  
No reutilizar  
Kertakäyttöinen  
Ne pas réutiliser  
Non riutilizzare  
Do jednorazowego użytku  
De uso único  
Återanvänds ej

**Literature** • Literatur • Bibliografia • Kirjallisuus • Bibliographie • Bibliografia • Bibliografia • Literature • Litteratur

**Heikkinen T.**

Comparative study of nasopharyngeal aspirate and nasal swab specimens for detection of influenza. BMJ (2001) 322:138.

**Ghebremedhin B. et al.**

Comparison of the performance of the rapid antigen detection actim Influenza A&B test and RT-PCR in different respiratory specimens. Journal of Medical Microbiology (2009), 58, 365-370.

**Steininger C. et al.**

Near-patient assays for diagnosis of influenza virus infection in adult patients. Journal Compilation 2009 European Society of Clinical Microbiology and Infectious Diseases, CMI, 15, 267-273.

**Stuchbery A. et al.**

Evaluation of a one step rapid immunochromatographic test for influenza A and B viruses. Australian Society for Microbiology National Conference in Canberra, Australia (25-29 September, 2005), Poster.

# Actim<sup>®</sup>

## Influenza A&B

Actim Oy  
Klovinpellontie 3, FI-02180 Espoo, Finland  
Tel. +358 9 547 680  
[www.medixbiochemica.com](http://www.medixbiochemica.com)  
[medix@medixbiochemica.com](mailto:medix@medixbiochemica.com)

© Actim Oy. All rights reserved.  
Actim is a trademark of Actim Oy.

AOACE32836-1 2019-04-25