
INSTRUCTIONS FOR USE


EN	Instructions for use	3
CS	Návod k použití	7
DA	Brugervejledning	11
DE	Gebrauchsanweisung	15
ES	Instrucciones de uso	19
FI	Käyttöohje	23
FR	Notice d'utilisation	27
IT	Istruzioni d'uso	31
NL	Gebruiksaanwijzing	35
NO	Bruksanvisning	39
PL	Instrukcja użytkowania	43
PT	Instruções de utilização	47
RO	Instrucțiuni de utilizare	51
RU	Инструкция по применению	55
SV	Bruksanvisning	59

Structure of dipstick

Test procedure and results

Structure of cartridge

Instructions for use

Numbers 1 - 13 refer to illustrations on inner cover.

Structure of dipstick

1 Dip area 2 Result area 3 Test line 4 Control line

Structure of cartridge

5 QR code

Intended use

The Actim® PROM 1ngeni test is a qualitative immunochromatographic dipstick test for the detection of amniotic fluid in vaginal secretions during pregnancy. Actim PROM 1ngeni detects IGFBP-1, which is a major protein in amniotic fluid and a marker of the amniotic fluid in a vaginal sample. The test is intended for professional use to help diagnose the rupture of fetal membranes (ROM) in pregnant women.

Kit components

Actim PROM 1ngeni kit (30831RETAC) contains 10 test packs with instructions for use. The components of each Actim PROM 1ngeni test pack (30821RETAC) are:

- One sterile polyester swab for specimen collection.
- One tube of Specimen Extraction Solution (0.5 ml). This phosphate-buffered solution contains bovine serum albumin (BSA), protease inhibitors and preservatives.
- One dipstick in a sealed aluminum foil pouch with desiccant.
- One cartridge for reading the result with the Actim 1ngeni instrument. The QR code on the cartridge includes information about the test (Actim PROM 1ngeni), the dipstick lot, expiry date, and lot specific calibration data.

Storage

Store the test kit at +2... +25° C. Stored unopened, each component can be used until the expiry date marked on the component. The kit can also be stored for 2 months at +2... +30° C. Use the dipsticks shortly after their removal from the aluminum foil pouch.

Specimen collection and sample handling

The specimen consists of vaginal secretion that is extracted into the Specimen Extraction Solution provided. A sample is obtained using a sterile polyester swab (provided in the test pack).

1 The sample should be collected prior to performing a digital examination and/or transvaginal ultrasound. 1 Take care not to touch anything with the swab before taking the sample.

1. Separate the labia and carefully insert the tip of the swab into the vagina toward the posterior fornix until resistance is met. Alternatively, the sample can be taken from the posterior fornix during a sterile speculum examination.
2. Leave the swab in the vagina for 10-15 seconds to allow it to absorb the vaginal secretion 6.
3. Open the Specimen Extraction Solution tube and put it in a vertical position. The specimen is extracted immediately from the swab by swirling the swab vigorously in the extraction solution for 10-15 seconds 7. Press the swab against the wall of the Specimen Extraction Solution tube to remove any remaining liquid from the swab. Discard the swab after extraction. **NOTE!** Do not leave the swab to the tube.

1 Specimens should be tested as soon as possible after extraction but, in any case, no more than 4 hours after the specimen collection and extraction. If a specimen cannot be tested within this time it should be frozen. After thawing, the specimens should be mixed and tested as described below.

Test procedure – Patient Test

In this procedure, test reaction occurs inside the instrument and results are read automatically at 5 minutes. Insertion of the cartridge into the Actim 1ngeni automatically initiates the Actim PROM 1ngeni test procedure.

1 If the test kit is stored refrigerated, allow the test pack to reach room temperature.

Prepare the instrument:

1. Switch the Actim 1ngeni on (if power off). After the automatic self test, log onto the instrument by entering your user ID either manually or with a barcode reader.
2. Select Test icon 8 and select Patient Test 9.

1 Do not choose [Patient test, Read only]. See below "Test procedure – read only".

3. Enter patient ID either manually or with a barcode reader.

Activate the test and start the procedure:

1. When the Actim 1ngeni screen shows "Activate and insert test", open the foil pouch containing the dipstick.
- 1 Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper turquoise part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch.

2. Activate the test by placing the yellow dip area of the test into the extracted sample **10** and hold it there until you see the liquid front enter the result area **11**.
 3. Immediately remove the dipstick from the solution, place it in the cartridge (reading area forward and coloured side facing up) and insert the cartridge inside the Actim 1ngeni (QR code forward and facing up) until the cartridge clicks into place **12**. Follow and respond to any information on the screen.
- !** Use only the cartridge included in the test pack. **!** Ensure that the dipstick is firmly inside the cartridge. **!** Ensure that the timer appears on the screen.

Proceed to section “Results”.

Test procedure – Patient Test, Read Only

In the read only-procedure, the test reaction occurs outside Actim 1ngeni. The user is responsible for the timing and for starting the analysis. Insertion of the cartridge into Actim 1ngeni automatically initiates the reading of the test result.

- !** As the test must be read at 5 minutes, ensure that Actim 1ngeni is available for analysis before activating the test. **!** If the test kit is stored in refrigeration, allow the test pack to reach room temperature.

Prepare the instrument:

1. Switch the Actim 1ngeni on (if power off). After the automatic self test, log onto the instrument by entering your user ID either manually or with a barcode reader.
2. Select **Test icon 8** and select **Patient Test, Read only**.
3. Enter the patient ID either manually or with a barcode reader.

Activate the test and start the procedure:

- !** Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper turquoise part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch.
4. When the Actim 1ngeni screen shows “Insert test”, open the foil pouch containing the dipstick.
5. Activate the test by placing the yellow dip area of the test into the extracted sample **10** and hold it there until you see the liquid front enter the result area **11**.
6. Immediately remove the dipstick from the solution, place it in the cartridge (reading area forward and coloured side facing up) and place the cartridge in a horizontal position and start your external timer for 5 minutes.

7. **Start the analysis:** Immediately at 5 minutes, insert the cartridge inside the Actim 1ngeni (QR code forward and facing up) until the cartridge clicks into place **12**. Follow and respond to any information on the screen until the “Analysing” screen appears.

- !** Use only the cartridge included in the test pack. **!** Ensure that the dipstick is firmly inside the cartridge.

Proceed to section “Results”.

Results

The Actim 1ngeni automatically reports and saves the result.

- !** Do not visually interpret the result under any circumstances after removing the test from the instrument because the test is optimised for the Actim 1ngeni.
- When the test result is visible on the screen, remove and discard the dipstick and the cartridge. The result can be either printed or exported into a USB stick or stored in the Actim 1ngeni for later use. Refer to the instrument manual for details.
 - Log off when leaving the instrument.

Positive and negative results with additional test-related information are shown on the screen **13**. If problems have occurred during the test process, an error is reported and no result is shown.

Quality control and Calibration

Test calibration and quality control

Each Actim PROM 1ngeni test kit lot has been calibrated and quality controlled during manufacturing to give reproducible results with the Actim 1ngeni instrument. QR code on the cartridge contains the lot specific calibration data.

Depending on local hospital practices and the local administrator for the Actim 1ngeni instrument, quality control for the Actim PROM 1ngeni test system using liquid controls may be advised. Refer to Actim PROM controls (30800ETAC) for details. There are no daily requirements for controlling the instrument performance using liquid controls.

Actim 1ngeni instrument

The Actim 1ngeni instrument has various levels of quality control procedures available. Refer to Actim 1ngeni instrument manual for details.

Limitations of the test

- The test is intended for *in vitro* diagnostic use only.
- If rupture of fetal membranes has occurred **but the leakage of amniotic fluid has ceased** more than 12 hours before the specimen is taken, IGFBP-1 may have been degraded by proteases in the vagina and the test may give a negative result.
- A negative test result is an indication on the present condition and cannot be used to predict the forthcoming.

Notes

- The test requires approximately 150 µl of extracted sample to ensure proper performance of the test.
- Use care when placing the dipstick in the sample tube. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing.
- Do not use the dipstick if its aluminum foil pouch or the seals of the pouch are not intact.
- Use only the swab provided with the Actim PROM kit.
- Improper sampling may lead to false result.
- When dipping, be careful to hold the dipstick in position (with the dip area in the sample extract) until the sample liquid front reaches the result area.
- The availability of the workflow and read only modes is dependent on the administrator settings.
- A test must always be read in read only mode in case the test reaction has occurred outside the instrument.
- Do not insert anything other than the cartridges provided with Actim 1ngeni tests into the instrument.
- Proper timing of reading the dipstick is essential. A delayed reading of the dipstick can have an effect on the test results. In these cases the test should be repeated.
- Ensure that the dipstick is firmly inside the cartridge and does not slip out when placed into the reader.
- If the procedural control line is not detected by the Actim 1ngeni instrument, the result will be declared indeterminate and the test should be repeated.

- The result must not be visually inspected after analysing it with the instrument.
- Do not re-analyse a test with the instrument.
- For additional information, please refer to the Actim 1ngeni instrument manual.
- A positive Actim PROM test result, although detecting the presence of amniotic fluid in the sample, does not locate the site of the rupture.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

Principle of the test

The concentration of IGFBP-1 (insulin-like growth factor binding protein-1) in amniotic fluid is 100 to 1,000 times higher than in maternal serum. IGFBP-1 is not usually present in the vagina, but after rupture of fetal membranes, amniotic fluid with a high concentration of IGFBP-1 mixes with vaginal secretions. In Actim PROM test, a specimen of vaginal secretion is taken with a sterile polyester swab and the specimen is extracted into Specimen Extraction Solution. The presence of IGFBP-1 in the solution is detected using a dipstick.

The test is based on immunochromatography. It uses two monoclonal antibodies to human IGFBP-1. One is bound to blue latex particles (the detecting label). The other is immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result. When the dip area of the dipstick is placed in an extracted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains IGFBP-1 it binds to the antibody labeled with latex particles. The particles are carried by the liquid flow and, if IGFBP-1 is bound to them, they bind to the catching antibody. Actim 1ngeni instrument will recognize a blue line (test line) if the concentration of IGFBP-1 in the sample exceeds the detection limit of the test. A second blue line detected by the instrument, the control line, confirms correct performance of the test.

Performance of the test

Analytical Sensitivity

The analytical sensitivity of Actim PROM 1ngeni was evaluated using samples with different concentrations of IGFBP-1. The samples were tested using three replicates. The analytical sensitivity (detection limit) is approximately 25 µg/l for Actim PROM 1ngeni test and the results remain positive at least to 500,000 µg/l.

Reproducibility

Reproducibility was assessed using samples containing 0 – 50 µg/L of IGFBP-1. Two operators performed the testing on different days. Three lots of the test were used; each lot read on a different instrument using the two test modes. The specimens were tested using ten replicates for a total of 120 tests per concentration. Negative samples were all identified negative and positive samples were all identified positive. One hundred percent (100 %) agreement with expected test results was obtained.

Method comparison

Performance of Actim PROM 1ngeni was compared to Actim PROM using spiked vaginal samples containing different concentration levels of IGFBP-1. The testing was performed by five operators using five different instruments. Overall agreement between the test formats was 100 %. These results confirm that the performance of Actim PROM 1ngeni test is similar to Actim PROM test (Table).

TABLE: METHOD COMPARISON BETWEEN ACTIM PROM 1NGENI AND ACTIM PROM

	Actim PROM POS	Actim PROM NEG	Total
Actim PROM 1ngeni POS	50	0	50
Actim PROM 1ngeni NEG	0	50	50
Total	50	50	100

Diagnostic performance

Actim PROM test has been evaluated in several clinical studies. Selection of typical studies is shown in FIG 1 on the inner back cover.

Interfering substances

Whole blood, semen, urine, pH, several drugs, shower and bath products, odor control products, and vaginal pathogens were tested with Actim PROM test and were found not to affect Actim PROM test performance.

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© Actim Oy. All rights reserved.
Actim is a trademark of Actim Oy.
Actim Oy