
INSTRUCTIONS FOR USE





Actim[®] Partus

EN	Instructions for use	3
CS	Návod k použití	6
DA	Brugervejledning	9
DE	Gebrauchsanweisung	12
ES	Instrucciones de uso	15
FI	Käyttöohje	18
FR	Notice d'utilisation	21
IT	Istruzioni d'uso	24
NL	Gebruiksaanwijsties	27
NO	Bruksanvisning	30
PL	Instrukcja użytkowania	33
PT	Instruções de utilização	36
RO	Instrucțiuni de utilizare	39
RU	Инструкция по применению	42
SV	Bruksanvisning	46

Structure of dipstick

Test procedure and results

Instructions for use

Numbers ① - ⑩ refer to illustrations on inner cover.

Structure of dipstick

① Dip area ② Result area ③ Test line ④ Control line

Intended use

Actim[®] Partus test is a visually interpreted, qualitative immunochromatographic dipstick test for detecting the presence of phosphorylated IGFBP-1 (insulin-like growth factor binding protein-1) in cervical secretions during pregnancy. The test is intended for professional use to help predict the risk of preterm or imminent delivery when fetal membranes are intact. A negative test result is a clear indication that the patient will not deliver within 7-14 days.

Kit components

Actim Partus kit 31931ETAC contains 10 test packs with instructions for use. The kit 31930ETAC contains one test pack and instructions for use. The components of each Actim Partus test pack (31921ETAC) are:

- One sterile polyester swab for specimen collection.
- One tube of Specimen Extraction Solution (0.5 ml). This phosphate-buffered solution contains bovine serum albumin (BSA), protease inhibitors and preservatives.
- One dipstick in a sealed aluminium foil pouch with desiccant.

Storage

Store the test kit at +2...+25 °C. Stored unopened, each component can be used until the expiry date marked on the component. The kit can also be stored for 2 months at +2...+30 °C. Use the dipsticks shortly after their removal from the aluminum foil pouch.

Specimen collection

The specimen is cervical secretion that is extracted into the Specimen Extraction Solution provided. The sample should be collected prior to performing digital examination and/or transvaginal ultrasound. A cervical secretion sample is obtained using a sterile polyester swab (provided in the kit) from the cervical os during a sterile speculum examination. Take care not to touch anything with the swab before taking the sample. The swab should be left in the cervical os for 10-15 seconds to allow it to absorb the secretion specimen ⑤.

Open the Specimen Extraction Solution tube and put it in a vertical position. The specimen is extracted immediately from the swab by swirling the swab vigorously in the extraction solution for 10-15 seconds ⑥. Press the swab against the wall of the Specimen Extraction Solution tube to remove any remaining liquid from the swab. Discard the swab after extraction. **NOTE!** Do not leave the swab to the tube.

Specimens should be tested as soon as possible after extraction but in any case no more than 4 hours after specimen collection and extraction. If a specimen cannot be tested within this time it should be frozen. After thawing, the specimens should be mixed and tested as described below.

Test procedure and interpretation of the results

1. If stored refrigerated, allow the aluminum foil pouch and the Extraction Solution tube to reach room temperature. Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper purple part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch.
2. Place the yellow dip area into the extracted sample ⑦ and hold it **there until you see the liquid front enter the result area ⑧**. Remove the dipstick from the solution and place it in a horizontal position.
3. The result can be interpreted as positive as soon as two blue lines become visible in the result area. Negative result should be read at **5 minutes ⑨**. **Do not pay attention to any lines appearing later than 5 minutes.**
4. If **two blue lines, the test line and control line**, appear, the test result is **positive**. If **one blue line, the control line**, appears, the test result is **negative**. **If the control line does not appear, the test is invalid ⑩.**

Limitations of the test

- The test is intended for *in vitro* diagnostic use only.
- Before performing the test ensure that the fetal membranes are intact (for example with Actim PROM test), because with ruptured fetal membranes Actim Partus test will also give a positive result.
- Patients with moderate or heavy vaginal bleeding should not be tested. It is recommended that a sample be taken when bleeding has stopped, and the extract is essentially blood-free.
- The test result indicates the risk at the time of sampling and changes in the patient's condition may later affect the final outcome of the pregnancy.

Notes

- No quantitative interpretation should be made based on the test results.
- The test requires about 150 µl of extracted sample to ensure proper performance of the test.
- Care must be taken when placing the dipstick in the sample tube. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue colouring in the result area before testing.
- Do not use the dipstick if its aluminium foil pouch or the seals of the pouch are not intact.
- Use only the swab provided with the kit.
- Improper sampling may lead to false negative result.
- When dipping, be careful to hold the dipstick in position (with the dip area in the sample extract) until the sample liquid front reaches the result area.
- If the control line does not appear, the test is invalid, and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly it is recommended that the test be repeated.

- At five minutes the appearance of any faint-to-dark blue test line along with a control line indicates a positive result. However, do not pay attention to any lines appearing after 5 minutes.
- If only the control line is visible, the result should be interpreted as negative only after 5 minutes have elapsed.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

Principle of the test

Decidual cells synthesize the phosphorylated forms of IGFBP-1 (pIGFBP-1) while amniotic fluid contains substantial quantities of non- and less phosphorylated forms of IGFBP-1. When delivery is approaching, fetal membranes begin to detach from the decidua parietalis, and small amounts of pIGFBP-1 begin to leak into cervical secretions. In Actim Partus test a cervical specimen sample is taken with a sterile polyester swab during sterile speculum examination and extracted into the Specimen Extraction Solution. The presence of pIGFBP-1 in the solution is detected using a dipstick.

The test is based on immunochromatography. It involves two monoclonal antibodies to human IGFBP-1. One is bound to blue latex particles (the detecting label). The other is immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result. When the dip area of the dipstick is placed in an extracted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains pIGFBP-1 it binds to the antibody labeled with latex particles. The particles are carried by the liquid flow and, if pIGFBP-1 is bound to them, they bind to the catching antibody. A blue line (test line) will appear in the result area if the concentration of pIGFBP-1 in the sample exceeds the detection limit of the test. A second blue line (control line) confirms correct performance of the test.

Performance of the test

Analytical Sensitivity

The analytical sensitivity (detection limit) of the Actim Partus test was identified by evaluating different concentrations of IGFBP-1 in extracted samples on three different lots of the Actim Partus test.

Two different operators each interpreted ten devices run at each concentration under various lighting conditions for a total of 60 determinations per level. The Actim Partus test limit of detection is approximately 10 µg/l in extracted sample. The measuring range of the Actim Partus is approximately 10-10 000 µg/l in extracted sample.

Analytical Specificity

Analytical specificity (cross-reactivity) was tested with human IGFBP proteins at a concentration level of 5000 µg/l of IGFBP-2, -4, -5 and -6 proteins, and at 50 000 µg/l of IGFBP-3 protein using one lot of the Actim Partus test. No cross-reactivity was seen using human IGFBP-2, -3, -4, -5 and -6 proteins. The Actim Partus test is specific to human IGFBP-1.

Repeatability

Two panels of specimens consisting of samples of different IGFBP-1 concentration levels were evaluated for intra-assay precision. Two different operators each interpreted one panel of specimens with 10 replicates during the same day using three different lots of the Actim Partus test for a total of 60 determinations per level. Repeatable results were obtained.

Diagnostic performance

Actim Partus test has been evaluated in several clinical studies. Selections of typical studies are shown in FIG 1 and in FIG 2 on the inner back cover. In these studies, Actim Partus test was used in predicting of pre-term labor within 7 days (FIG 1) and in predicting of labor before 32-37 weeks (FIG 2).

Interference testing

The following substances, conditions and micro-organisms were tested with the Actim Partus test and were found not to affect the Actim Partus test performance.

TESTED SUBSTANCE/CONDITION/ MICRO-ORGANISM	INTERFERENCE (NO/YES)
Semen	No
Urine	No
Surgical lubricant	No
Betadine	No
Personal Lubricant	No
Hexochlorophene	No
Hexochlorophene cream	No
Econazole cream	No
Prostaglandin E2	No
Ampicillin	No
Erythromycin	No
Dexamethasone	No
Betamethasone	No
Naproxen	No
Nifedipine	No
Magnesium sulfate	No
Oxytocin	No
Terbutaline	No
17-alpha hydroxyprogesterone caproate	No
Progesterone capsule (oral use)	No
Progesterone capsule (vaginal use)	No
Indomethacin	No
Bilirubin	No
Shower and bath products	No
Baby Oil	No
Baby Powder	No
Feminine Deodorant Suppositories	No
Vaginal Gel	No
Feminine Deodorant Film	No
Candida albicans	No
Gardnerella vaginalis	No
Neisseria gonorrhoea	No
Chlamydia trachomatis	No
HSV-1	No
HSV-2	No

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Akercan F et al.

Value of cervical phosphorylated insulinlike growth factor binding protein-1 in the prediction of preterm labor. J Reprod Med (2004) 49: 368-372.

Altinkaya O et al.

Cervical phosphorylated insulin-like growth factor binding protein-1 in prediction of preterm delivery. Arch Gynecol Obstet. (2009) 279:279-283.

Azlin MI et al.

Role of pIGFBP-1 and ultrasound cervical length in predicting pre-term labour. Journal of Obstetrics and Gynaecology (2010) 30:456-460.

Brik Spinelli M et al.

Phosphorylated insulin-like growth factor binding protein-1 and cervical measurement in women with threatening preterm birth. Acta Obstet Gynecol Scand (2010) 89:268-74.

Eroglu D et al.

Prediction of preterm delivery among women with threatened preterm labor. Gynecol Obstet Invest (2007) 64:109-116.

Kekki M et al.

Insulin-like growth factor-binding protein-1 in cervical secretion as a predictor of preterm delivery. Acta Obstet Gynecol Scand (2001) 80:546-551.

Lembet A et al.

New rapid bed-side test to predict preterm delivery: phosphorylated insulin-like growth factor binding protein-1 in cervical secretions. Acta Obstet Gynecol Scand (2002) 81:706-712.

Rutanen EM

Insulin-like growth factors in obstetrics. Opin Obstet Gynecol (2000) 12:163-168.

Tanir HM , Sener T, Yildiz Z.

Cervical phosphorylated insulin-like growth factor binding protein-1 for the prediction of preterm delivery in symptomatic cases with intact membranes. J Obstet Gynaecol Res (2009) 1:66-72.

Ting HS et al.

Comparison of bedside test kits for prediction of preterm delivery: phosphorylated insulin-like growth factor binding protein-1 (pIGFBP-1) test and fetal fibronectin test. Ann Acad Med Singapore (2007) 36:399-402.

Actim[®] Partus

Actim Oy
Klovinpellontie 3, FI-02180 Espoo, Finland
Tel. +358 9 547 680
www.medixbiochemica.com
medix@medixbiochemica.com

© Actim Oy. All rights reserved.
Actim is a trademark of Actim Oy.

AOACE31931-1 2019-04-24