

THE QUANTITATIVE ACTIM® PARTUS 1NGENI RAPID TEST RELIABLY PREDICTS PRETERM OR IMMINENT DELIVERY IN SYMPTOMATIC PREGNANT WOMEN

H. Miller¹, H. How², R. Newman³, M. Nageotte⁴, M. Fries⁵, R. Ramsey⁶, C. Zavala⁶, L. Kanto⁷, E. Svens⁷, J. Juhila⁷

1. Watching Over Mother and Babies Foundation, 5301 E. Grant Road, Tucson, Arizona 85712, USA
2. Norton Hospital "Children Foundation Building, 601 South Floyd Street, Louisville, Kentucky 40202, USA
3. Medical University of South Carolina, 96 Jonathan Lucas Street, Charleston, South Carolina 29425, USA
4. Long Beach Memorial Medical Center, 2801 Atlantic Ave, Long Beach, California 92801, USA
5. Washington Hospital Center, Department of Obstetrics and Gynecology, 110 Irving Street NW, 5B-11, Washington, DC 20010, USA
6. Methodist University Hospital, 1265 Union Avenue 8 East, Memphis, Tennessee 38104, USA
7. Actim Oy, Klövinpellontie 3, FI-02180 Espoo, Finland

Elevated cervical pHIGFBP-1 indicates preterm and imminent birth

Spontaneous preterm birth (sPTB) continues to be responsible for the majority of neonatal morbidity and mortality. Only about 20% of women presenting with suspected preterm labour will actually deliver preterm.¹ In order to be able to focus treatment, it is important to have a test to help predict who is most likely to deliver preterm.

An elevated concentration of the highly phosphorylated isoform of insulin-like growth factor binding protein 1 (pHIGFBP-1), produced in decidua^{2,3,4}, in cervical secretions is a predictor of preterm or imminent birth.^{5,6}

Actim® Partus 1ngeni is an immunochromatographic, quantitative rapid test that employs monoclonal antibodies specific for decidual pHIGFBP-1 to identify women with intact membranes who are at increased risk for preterm birth.^{5,7,8,9,10} We retrospectively evaluated the clinical performance of the quantitative Actim Partus 1ngeni to predict sPTB within 7 and 14 days of sample collection at different gestational ages.

Materials & Methods

Cervical pHIGFBP-1 concentration from 691 samples obtained at hospitals or maternal clinics was quantified using the fully automated Actim® 1ngeni System (Figure 1) with quantitative Actim Partus 1ngeni test according to manufacturer's instructions. The key inclusion criteria were: singleton pregnancy with signs and/or symptoms of labor, gestational age ≥ 22 and < 36 weeks, and intact amniotic membranes. Women with heavy vaginal bleeding were excluded from the study. The primary clinical outcome of the study was the occurrence of sPTB after sample collection. The relative risk (RR) of sPTB within 7 and 14 days was calculated for elevated pHIGFBP-1 concentration categories defined in the Actim Partus 1ngeni test, as compared to pHIGFBP-1 concentrations of ≤ 2 $\mu\text{g/mL}$.

Figure 1. pHIGFBP-1 concentration was quantified from cervical swab samples using an automated and digital Actim 1ngeni System with Actim Partus 1ngeni test.



Elevated Actim Partus 1ngeni test result reliably predicts delivery within 7 and 14 days

Of the women included in this study, 360 (52%) experienced a spontaneous labor. Nine of these women (3%) delivered within 7 days and 21 (6%) within 14 days of sample collection. Among women with cervical pHIGFBP-1 concentrations of 50.0–249.9 $\mu\text{g/mL}$ and ≥ 250 $\mu\text{g/mL}$ the RR of sPTB within 7 days of sample collection was 15.8 ($p < 0.05$) and 22.8 ($p < 0.01$), respectively. Additionally, the RR of sPTB within 14 days of sample collection was increased for women with cervical pHIGFBP-1 concentrations of 50.0–249.9 $\mu\text{g/mL}$ (4.5; $p < 0.05$) and ≥ 250 $\mu\text{g/mL}$ (8.1; $p < 0.001$).

Table 1. The increasing concentration of cervical pHIGFBP-1 correlates with increased relative risk (RR) of delivery within 7 and 14 days from sample collection.

pHIGFBP-1 (ng/mL)	Delivery ≤ 7 days	Delivery ≤ 14 days
≤ 2	-	-
2.1-9.9	5.9	0.8
10.0-49.9	10.0	2.1
50.0-249.9	15.8*	4.5*
≥ 250	22.8**	8.1***

Relative risk compared to pHIGFBP-1 ≤ 2 ng/mL
z-statistics significance: * $P < 0.05$, ** $P < 0.01$, *** $P < 0.001$

Conclusions

In this quantitative retrospective analysis, quantitative Actim Partus 1ngeni results with increasing cervical pHIGFBP-1 concentrations are associated with increasing risk of sPTB within 7 and 14 days of sample collection. The Actim Partus 1ngeni holds promise as a rapid point-of-care test for identifying women at increased risk of sPTB. This detection of different risk groups allows for more directed interventions for high-risk patients and reduces the cost and unnecessary treatment for patients in low risk of sPTB.

Acknowledgements: We wish to thank our Laboratory Technicians for their excellent technical assistance in this study.

References

1. Goldenberg et al. The preterm prediction study: fetal fibronectin testing and spontaneous preterm birth. NICHD Maternal Fetal Medicine Units Network. *Obstet Gynecol* 1996;87:643-8.
2. Rutanen et al. Radioimmunoassay of placental protein 12: Levels in amniotic fluid, cord blood, and serum of healthy adults, pregnant women and patients with trophoblastic disease. *Am J Obstet Gynecol* 1982;144:460-63.
3. Julkunen et al. Primary structure of human insulin-like growth factor binding protein/placental protein 12 and tissue expression of its mRNA. *FEBS Lett* 1988;236:295-02.
4. Martina et al. Gestational age-dependent expression of insulin-like growth factor binding protein-1 (IGFBP-1) phosphoisoforms in human extra-amniotic cavities, maternal serum, and decidua suggests decidua as the primary source of IGFBP-1 in these fluids during early pregnancy. *J Clin Endocrinol Metab* 1997;82:1894-98.
5. Kekki et al. Insulin-like growth factor-binding protein-1 in cervical secretion as a predictor of preterm delivery. *Acta Obstet Gynecol Scand* 2001;80:546-51.
6. Paternoster et al. Cervical pHIGFBP-1 in the evaluation of the risk of preterm delivery. *Acta Obstetrica et Gynecologica* 2007;86:151-55.
7. Nuutila et al. Phosphorylated isoforms of insulin-like growth factor binding protein-1 in the cervix as a predictor of cervical ripeness. *Obstetrics & Gynecology*, 1999; 94(2), 243-249.
8. Tanir et al. Cervical phosphorylated insulin-like growth factor binding protein-1 for the prediction of preterm delivery in symptomatic cases with intact membranes. *J Obstet Gynaecol Res* 2009;35(1):66-72.
9. Bruijn et al. Comparison of the Actim Partus test and the fetal fibronectin test in the prediction of spontaneous preterm birth in symptomatic women undergoing cervical length measurement. *Eur J Obstet Gynecol Reprod Biol* 2016;206:220-224.
10. Vallikkannu et al. Insulin-like growth factor binding protein 1, Bishop score, and sonographic cervical length: tolerability and prediction of vaginal birth and vaginal birth within 24 hours following labour induction in nulliparous women. *BJOG* 2016; DOI: 10.1111/1471-0528.14175.