INSTRUCTIONS FOR USE

Actim[®] SARS-CoV-2

EN	Instructions for use	3
CS	Návod k použití	6
DA	Brugervejledning	9
DE	Gebrauchsanweisung	12
ES	Instrucciones de uso	15
FI	Käyttöohje	18
NL	Gebruiksaanwijzing	21
NO	Bruksanvisning	24
PL	Instrukcja użytkowania	27
SV	Bruksanvisning	30

Structure of dipstick

[1] Dip area





[2] Result area

Test line [3]

[4] Control line

Test procedure and results

Instructions for use

Numbers [1] - [9] refer to illustrations on inner cover.

Structure of dipstick

[1] Dip area [2] Result area [3] Test line [4] Control line

Intended use

The Actim® SARS-CoV-2 is a visually interpreted, *in vitro* diagnostic immunochromatographic rapid test for the qualitative detection of SARS-CoV-2 virus antigens in nasopharyngeal swab samples. The test is intended for professional use to help diagnose acute COVID-19 infection.

Kit components

The components of the Actim SARS-CoV-2 kit 35522ETAC are:

- 20 individually packed dipsticks (35523ETAC) in sealed aluminum foil pouches with desiccant.
- 20 Swab Extraction Buffer tubes (35524ETAC) needed for sample extraction. Each tube contains 0.3ml of buffer, which contains detergents, bovine serum albumin (BSA) and preservatives.
- 20 swabs for swab sample collection.
- Instructions for use

Materials required but not provided: Clock/Timer, Personal protective equipment

Storage

Store the test kit at +2...+30 °C. Stored unopened, each component can be used until the expiry date marked on the component. Use the dipsticks shortly after their removal from the aluminum foil pouch.

Specimen collection and handling

Collect a nasopharyngeal swab sample [5]. Gently insert the swab horizontally into the nostril of the patient, reaching the surface of posterior nasopharynx that present the most secretion under visual inspection. Keep the swab in place and rotate slowly five [5] times against the nasal wall and gently remove from nose. This should take approximately 10-15 seconds.

The samples should be tested as soon as possible after collection but in any case, no more than 1 hour after specimen collection and extraction.

Open the Swab Extraction Buffer tube and put it in a vertical position. The specimen is extracted immediately from the swab by swirling the swab vigorously in the extraction solution for 10–15 seconds [6]. Press the swab against the wall of the Swab Extraction Buffer tube to remove any remaining liquid from the swab. Discard the swab after extraction.

Note! Do not leave the swab to the tube.

Test procedure and interpretation of the results

1. Allow the aluminum foil pouch to reach room temperature if stored refrigerated. Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper blue part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch.
2. Place the yellow dip area [the lower end of the dipstick] into the test tube with extracted sample and leave it there until reading [7].
3. Read the result at 15 minutes. A result can be interpreted as positive as soon as one blue and one black line become visible in the result area. If only one black line appears within 15 minutes, the result is interpreted as negative [8]. Appearance of a black control line confirms that the test has been performed properly.

If, in addition to the black control line you see:

- **blue line** in the result area, the result is **SARS-CoV-2 positive**.
- **no blue line** in the result area, the result is **SARS-CoV-2 negative**.

If the black control line does not appear, the test is **invalid [9]**.

Note! The appearance of any faint-to-dark blue test line along with a black control line indicates a positive result.

Quality Control procedure

External positive control is sold separately, REF VP-S-19-02, Microbix.

Perform according to the test procedure.

Limitations of the test

- The test is intended for *in vitro* diagnostic use only.
- Failure to follow instructions for test procedure, using wrong specimen or wrong interpretation of test results may affect test performance or/and produce invalid results.
- Positive test result does not rule out co-infections with other pathogens.
- A confirmed diagnosis should only be made by a physician with all clinical data available.

Notes

- Test is for single use.
- The test requires approximately 150 µl of extracted sample to ensure proper performance of the test.
- No quantitative interpretation should be made based on the test results.
- Use only the swab and Swab Extraction Buffer provided with the kit.
- A negative test result may occur if the sample was collected or stored improperly.
- The test detects both viable and non-viable SARS-CoV-2 virus.
- Do not use the dipstick or swab if their foil pouches or seals of the pouches are not intact.
- Care must be taken when placing the dipstick in the sample tube. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue or black coloring in the result area before testing.
- If a control line does not appear, the test is invalid and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly it is recommended that the test be repeated with another dipstick.
- All biological specimens and materials must be treated as potentially hazardous and disposed of in accordance with local authority guidelines.

- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.

Principle of the test

The test includes the extraction and detection of SARS-CoV-2 virus antigens. The test is based on immunochromatography. It utilizes monoclonal antibodies to SARS-CoV-2 virus antigens. Antibodies against SARS-CoV-2 virus antigens are bound to blue latex particles (the detecting label). Antibodies against SARS-CoV-2 virus antigens are immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result.

When the dip area of the dipstick is placed in an extracted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains SARS-CoV-2 virus, the virus antigen binds to the antibodies labeled with latex particles. The particles are carried by the liquid flow and, if SARS-CoV-2 virus antigen is bound to them, they bind to the corresponding catching antibody on the membrane.

A blue line (test line) will appear in the result area if the concentration of SARS-CoV-2 virus in the sample exceeds the detection limit of the test.

Performance of the test

Analytical Sensitivity

The analytical sensitivity of Actim SARS-CoV-2 test was assessed using SARS-CoV-2 viral strain FIN-25. The detection limit of SARS-CoV-2 virus is 0.5×10^2 TCID₅₀/ml and the result remain positive at least to 2×10^5 TCID₅₀/ml.

Analytical Specificity

Actim SARS-CoV-2 was evaluated with the viruses listed below. A negative test result was obtained with all viruses. Test results indicate that Actim SARS-CoV-2 test is specific for COVID-19.

VIRUS:	CONCENTRATION TESTED:
Human coronavirus 229E	5×10^4 TCID ₅₀ /ml
Human coronavirus OC43	1×10^9 TCID ₅₀ /ml
Influenza A virus	5×10^4 pfu/ml
Influenza B virus	2×10^4 pfu/ml
MERS	60 mg/l
Enterovirus	1×10^4 pfu/ml
Respiratory syncytial virus (A&B)	2×10^4 pfu/ml

Repeatability and Reproducibility

A panel of samples containing SARS-CoV-2 antigen (low positive and high positive) and a negative sample was evaluated for intra- (repeatability) and inter- (reproducibility) assay precision. For repeatability the samples were tested with 10 replicates during the same day using three lots of the Actim SARS-CoV-2 test. For reproducibility, the samples were tested by three operators on five days using three lots of the Actim SARS-CoV-2 test. Test shows consistent performance within run, between runs, between operators, between lots and between days. One hundred percent (100%) agreement with expected test results was obtained.

Interfering substances

Different nasal sprays [5%], nasal corticosteroid [5%] and whole blood [5%] were tested with Actim SARS-CoV-2 test. No interference of these substances or whole blood was observed with the test performance.

Diagnostic performance

Actim SARS-CoV-2 test was evaluated at University of Turku, Department of Virology, Finland. The results are shown in FIG 1. and FIG 2. on the inner back cover.

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© Actim Oy. All rights reserved.
Actim is a trademark of Actim Oy.
Actim Oy