INSTRUCTIONS FOR USE

Actim®
SARS-CoV-2

Actim®
Influenza A&B


EN	Instructions for use	3
DA	Brugervejledning	7
DE	Gebrauchsanweisung	11
ES	Instrucciones de uso	15
FI	Käyttöohje	19
FR	Notice d'utilisation	23
NO	Bruksanvisning	28
PL	Instrukcja użytkowania	32
SV	Bruksanvisning	36

Structure of dipstick

[1] Dip area

[2] Result area

Actim SARS-CoV-2

[5] Dip area








[6] Result area

Actim Influenza A&B

Test procedure and results

Instructions for use

Numbers [1] - [13] refer to illustrations on inner cover.

Structure of dipstick

Actim SARS-CoV-2

[1] Dip area [2] Result area [3] Test line [4] Control line

Actim Influenza A&B

[5] Dip area [6] Result area [7] Test lines [8] Control line

Intended use

Actim SARS-CoV-2

The Actim SARS-CoV-2 is a visually interpreted, *in vitro* diagnostic immunochromatographic rapid test for the qualitative detection of SARS-CoV-2 virus antigens in nasal or nasopharyngeal swab samples. The test is intended for professional use to help diagnose acute COVID-19 infection.

Actim Influenza A&B

The Actim Influenza A&B test is a visually interpreted, qualitative immunochromatographic rapid test for detection of influenza type A and type B virus antigens in respiratory samples. The test is intended for professional use to help diagnose acute influenza infections.

Kit components

The components of the Actim SARS-CoV-2 + Actim Influenza A&B kit 36031ETAC are:

- 10 individually packed dipsticks of Actim SARS-CoV-2 [35523ETAC] and Actim Influenza A&B [32801ETAC] in sealed aluminum foil pouches with desiccant.
- 10 Swab Extraction Buffer tubes [35524ETAC] needed for sample extraction. Each tube contains 0.3ml of buffer, which contains detergents, bovine serum albumin (BSA) and preservatives.
- 10 swabs for swab sample collection.
- Instructions for use.

Materials required but not provided: Clock/Timer, Personal protective equipment.

Storage

Store the test kit at +2...+30 °C. Stored unopened, each component can be used until the expiry date marked on the component. Use the dipsticks shortly after their removal from the aluminum foil pouch.

Specimen collection and handling

Collect a nasal [9a] or nasopharyngeal [9b] swab sample. Keep the swab in place and rotate slowly five (5) times against the nasal wall and gently remove from nose. This should take approximately 10-15 seconds.

The samples should be tested as soon as possible after collection but in any case, no more than 1 hour after specimen collection and extraction.

Open the Swab Extraction Buffer tube and put it in a vertical position. The specimen is extracted immediately from the swab by swirling the swab vigorously in the extraction solution for 10–15 seconds [10]. Press the swab against the wall of the Swab Extraction Buffer tube to remove any remaining liquid from the swab. Discard the swab after extraction.

Note! Do not leave the swab to the tube.

Test procedure

1. Allow the aluminum foil pouches to reach room temperature if stored refrigerated. Open the foil pouches containing the dipsticks by tearing. Do not touch the yellow dip area at the lower part of the dipsticks. Identifying marks may be written on the upper light and dark blue part of the dipsticks. The dipsticks must be used shortly after its removal from the foil pouch.
2. Place the yellow dip areas (the lower end of the dipsticks) backsides facing each other into the test tube with extracted sample, and leave those there until reading [11].

Interpretation of the result

The test lines are read at two different timepoints: at 10 minutes (Actim Influenza A&B) and at 15 minutes (Actim SARS-CoV-2).

Actim Influenza A&B (light blue dipstick)

Read the result at 10 minutes. A result can be interpreted as positive as soon as one or two blue and one black line become visible in the result area. If only one black line appears within 10 minutes, the result is interpreted as negative [12a]. Appearance of a black control line confirms that the test has been performed properly.

If, in addition to the black control line you see:

- **one blue line** in the **lower part** of the result area, the result is **influenza A positive**.
- **one blue line** in the **middle part** of the result area, the result is **influenza B positive**.
- **two blue lines** in the result area, the result is **influenza A and B positive**.
- **no blue lines** in the result area, the result is **influenza A and B negative**.

If the control line does not appear the test is invalid [12b].

Note! The appearance of any faint-to-dark blue test line along with a black control line indicates a positive result.

Actim SARS-CoV-2 (dark blue dipstick)

Read the result at **15 minutes**. A result can be interpreted as positive as soon as one blue and one black line become visible in the result area. If only one black line appears within 15 minutes, the result is interpreted as negative **[13a]**. Appearance of a black control line confirms that the test has been performed properly.

If, in addition to the black control line you see:

- **blue line** in the result area, the result is **SARS-CoV-2 positive**.
- **no blue line** in the result area, the result is **SARS-CoV-2 negative**.

If the black control line does not appear, the test is invalid [13b].

Note! The appearance of any faint-to-dark blue test line along with a black control line indicates a positive result.

Quality Control procedure

External positive controls are available separately.

Limitations of the test

- The test is intended for *in vitro* diagnostic use only.
- Failure to follow instructions for test procedure, using wrong specimen or wrong interpretation of test results may affect test performance or/and produce invalid results.
- Positive test result(s) does not rule out co-infections with other pathogens.

- A confirmed diagnosis should only be made by a physician with all clinical data available.

Notes

- Test is for single use.
- One Swab Extraction Buffer tube contains the volume sufficient for one Actim SARS-CoV-2 and one Actim Influenza A&B dipstick.
- No quantitative interpretation should be made based on the test results.
- Use only the swab and Swab Extraction Buffer provided with the kit.
- A negative test result may occur if the sample was collected or stored improperly.
- The test detects both viable and non-viable SARS-CoV-2 and influenza viruses.
- Do not use the dipstick or swab if their foil pouches or seals of the pouches are not intact.
- Care must be taken when placing the dipstick in the sample tube. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue or black coloring in the result area before testing.
- If a control line does not appear, the test is invalid and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly it is recommended to perform a new test.
- All biological specimens and materials must be treated as potentially hazardous and disposed of in accordance with local authority guidelines.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- The performance data of Actim Influenza is obtained with the Actim Influenza A&B test version 32836ETAC and thus contains data from all respiratory samples. If you wish to test influenza from other than nasopharyngeal sample, e.g. nasal swab sample, nasal aspirate and nasal wash sample, please order the kit 32836ETAC.

Principle of the test

Actim SARS-CoV-2

The test includes the extraction and detection of SARS-CoV-2 virus antigens. The test is based on immunochromatography. It utilizes monoclonal antibodies to SARS-CoV-2 virus antigens. Antibodies against SARS-CoV-2 virus antigens are bound to blue latex particles (the detecting label). Antibodies against SARS-CoV-2 virus antigens are immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result.

When the dip area of the dipstick is placed in an extracted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains SARS-CoV-2 virus, the virus antigen binds to the antibodies labeled with latex particles. The particles are carried by the liquid flow and, if SARS-CoV-2 virus antigen is bound to them, they bind to the corresponding catching antibody on the membrane.

A blue line (test line) will appear in the result area if the concentration of SARS-CoV-2 virus in the sample exceeds the detection limit of the test.

Actim Influenza A&B

The test includes the extraction and detection of influenza type A and B virus antigens. The test is based on immunochromatography. It utilizes monoclonal antibodies to influenza type A and type B virus antigens. Antibodies against both types are bound to blue latex particles (the detecting label). Antibodies against both types are immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result.

When the dip area of the dipstick is placed in an extracted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains influenza type A or B, the virus antigen binds to either of the antibodies labeled with latex particles. The particles are carried by the liquid flow and, if influenza virus antigen is bound to them, they bind to the corresponding catching antibody on the membrane. A blue line (test line) will appear in the lower part of the result area if the concentration of influenza type A in the sample exceeds the detection limit of the test. A blue line (test line) will appear in the middle of the result area if the concentration of influenza type B in the sample exceeds the detection limit of the test. Appearance of the black control line confirms correct performance of the test.

Performance of the test

Actim SARS-CoV-2

Analytical Sensitivity

The analytical sensitivity of Actim SARS-CoV-2 test was assessed using SARS-CoV-2 viral strain FIN-25. The detection limit of SARS-CoV-2 virus is 0.5×10^2 TCID₅₀/ml in extracted sample and the result remain positive at least to 2×10^5 TCID₅₀/ml.

Analytical Specificity

Actim SARS-CoV-2 was evaluated with the viruses listed below. A negative test result was obtained with all viruses. Test results indicate that Actim SARS-CoV-2 test is specific for COVID-19.

Viruses

	Concentration tested
Human coronavirus 229E	5×10^6 TCID ₅₀ /ml
Human coronavirus OC43	1×10^9 TCID ₅₀ /ml
Human coronavirus NL63	14.2 µg/ml
Influenza A virus	5×10^6 pfu/ml
Influenza B virus	2×10^6 pfu/ml
MERS	60 mg/l
Enterovirus	1×10^6 pfu/ml
Respiratory syncytial virus (A&B)	2×10^6 pfu/m
Respiratory syncytial virus	450 mg/l

Repeatability and Reproducibility

A panel of samples containing SARS-CoV-2 antigen (low positive and high positive) and a negative sample was evaluated for intra- (repeatability) and inter- (reproducibility) assay precision. For repeatability the samples were tested with 10 replicates during the same day using three lots of the Actim SARS-CoV-2 test. For reproducibility, the samples were tested by three operators on five days using three lots of the Actim SARS-CoV-2 test. Test shows consistent performance within run, between runs, between operators, between lots and between days. One hundred percent (100%) agreement with expected test results was obtained.

Interfering substances

Different nasal sprays (5%), nasal corticosteroid (5%) and whole blood (5%) were tested with Actim SARS-CoV-2 test. No interference of these substances or whole blood was observed with the test performance.

Diagnostic performance

Actim SARS-CoV-2 test was evaluated at University of Turku, Department of Virology, Finland. The results are shown in FIG 1. and FIG 2. on the inner back cover.

Actim Influenza A&B

Analytical Sensitivity

The analytical sensitivity of Actim Influenza A&B test was assessed using different influenza viral strains. The detection limit of influenza virus is 1.0×10^4 TCID₅₀/ml for influenza A/New Caledonia/20/99 [H1N1], 1.0×10^4 TCID₅₀/ml for influenza A/Panama/2007/99 [H3N2] and 2.5×10^4 TCID₅₀/ml for influenza B/Yamanashi/166/98. The test has further been shown to detect cultured avian influenza A subtypes [H5N1, H5N3, H7N3, H7N9, H9N2].

Analytical Specificity

Actim Influenza A&B test was evaluated with the microorganisms listed below. A negative test result was obtained with all microorganisms. Test results indicate that Actim Influenza A&B test is specific for influenza A and influenza B.

Viruses/Bacteria ($\geq 10,000$ TCID₅₀/ml)

- Echovirus types 2, 3, 6, 9, 11, 25, 30
- Adenovirus types 1, 2, 3, 4, 5, 6, 7, 10, 18
- Parainfluenza types 1, 2, 3
- Coxsackievirus types A9, B5, B6
- Rhinovirus types 1A, 2, 13, 15, 37
- Respiratory syncytial virus [RSV] types A and B
- Cytomegalovirus AD169
- Herpes simplex virus types 1 and 2
- Mycoplasma pneumonia ($>10,000$ CFU)
- SARS-CoV-2 nucleoprotein (10 µg/ml)
- Human coronavirus NL63 (14.2 µg/ml)

Bacteria/Yeast ($\geq 10^5$ cfu/ml)

- *Bordetella pertussis*
- *Bacteroides fragilis*
- *Neisseria gonorrhoeae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Candida albicans*
- *Streptococcus viridans*
- *Streptococcus pneumoniae*
- *Listeria monocytogenes*
- *Streptococcus pyogenes*
- *Streptococcus agalactiae*
- *Streptococcus* ssp. Group C
- *Streptococcus* ssp. Group G
- *Moraxella catarrhalis*
- *Enterococcus faecalis*
- *Staphylococcus Aureus ATCC 25923*
- *Staphylococcus epidermidis*
- *Klebsiella pneumoniae*
- *Pseudomonas aeruginosa*
- *Escherichia coli*
- *Proteus mirabilis*
- *Staphylococcus aureus ATCC 29213*
- *Staphylococcus aureus ATCC 43300*

Repeatability and Reproducibility

A panel of specimens consisting of samples of influenza A (low positive and high positive), samples of influenza B (low positive and high positive) and a negative sample was evaluated for intra- (repeatability) and inter- (reproducibility) assay precision. For repeatability the samples were tested with 10 replicates during the same day using one lot of the Actim Influenza A&B test. For reproducibility, the samples were tested on 7 days using one lot of the Actim Influenza A&B test. One hundred percent (100%) agreement with expected test results was obtained.

Interfering Substances

Nasal drops (30 %) and ibuprofen (15 mg/ml) were tested with Actim Influenza A&B test. No interference of these substances was observed with the test performance. Whole blood was tested with Actim Influenza A&B test. No interference of whole blood (2%) was observed with the test performance.

Diagnostic Performance

Actim Influenza A&B test was evaluated at Institute of Medical & Veterinary Science, Adelaide, Australia and at University of Turku, Department of Virology, Finland. The results are shown in FIG 3. and FIG 4. on the inner back cover.

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- **Heikkinen T.** Comparative study of nasopharyngeal aspirate and nasal swab specimens for detection of influenza. BMJ [2001] 322:138.
- **Ghebremedhin B. et al.** Comparison of the performance of the rapid antigen detection actim Influenza A&B test and RT-PCR in different respiratory specimens. Journal of Medical Microbiology [2009], 58, 365-370.
- **Steininger C. et al.** Near-patient assays for diagnosis of influenza virus infection in adult patients. Journal Compilation 2009 European Society of Clinical Microbiology and Infectious Diseases, CMI, 15, 267-273.
- **Stuchbery A. et al.** Evaluation of a one step rapid immunochromatographic test for influenza A and B viruses. Australian Society for Microbiology National Conference in Canberra, Australia [25-29 September, 2005], Poster.

FIG 1.

RESULTS FROM NASOPHARYGEAL SWAB SAMPLE

Actim SARS-CoV-2 Reference method: RT-PCR

	POS	NEG	TOTAL
POS	85	0	85
NEG	15	149	164
TOTAL	100	149	249

RESULTS FROM NASAL SWAB SAMPLE

Actim SARS-CoV-2 Reference method: RT-PCR

	POS	NEG	TOTAL
POS	32	0	32
NEG	4	44	48
TOTAL	36	44	80

FIG 2.

RESULTS FROM NASOPHARYGEAL SWAB SAMPLE

	Ct < 25	Ct < 29	Ct < 33	All samples
Sensitivity	100%	100%	93%	85%
Specificity				100%
True positives	44/44	63/63	78/84	85/100

RESULTS FROM NASAL SWAB SAMPLE

	Ct < 25	Ct < 29	Ct < 33	All samples
Sensitivity	100%	100%	94%	89%
Specificity				100%
True positives	21/21	26/26	32/34	32/36

FIG 3.

INFLUENZA A			REFERENCE METHOD EIA			INFLUENZA B			REFERENCE METHOD EIA		
Actim Influenza A&B	POS	NEG	TOTAL	Actim Influenza A&B	POS	NEG	TOTAL	Actim Influenza A&B	POS	NEG	TOTAL
POS	16	3	19	POS	0	1	1	POS	0	1	1
NEG	0	269	269	NEG	0	260	260	NEG	0	260	260
TOTAL	16	272	288	TOTAL	0	261	261	TOTAL	0	261	261
95% CI			95% CI			95% CI			95% CI		
SENS	100%	76–100%		SENS*				SPEC	100%	98–100%	
SPEC	99%	97–100%		SPEC				PPV	84%	60–96%	
PPV	84%	60–96%		PPV*				NPV	100%	98–100%	
NPV	100%	98–100%		NPV*							

*Due to limited number of positive influenza B infections in the study population, it is impossible to define the clinical performance.

FIG 4.

INFLUENZA A			REFERENCE METHOD TR-FIA			INFLUENZA B			REFERENCE METHOD TR-FIA		
Actim Influenza A&B	POS	NEG	TOTAL	Actim Influenza A&B	POS	NEG	TOTAL	Actim Influenza A&B	POS	NEG	TOTAL
POS	18	0	18	POS	19	0	19	POS	19	0	19
NEG	3	68	71	NEG	2	68	70	NEG	2	68	70
TOTAL	21	68	89	TOTAL	21	68	89	TOTAL	21	68	89
95 % CI			95 % CI			95 % CI			95 % CI		
SENS	86%	63–96%		SENS				SPEC	90%	68–98%	
SPEC	100%	93–100%		SPEC	100%	93–100%		PPV	100%	79–100%	
PPV	100%	78–100%		PPV	100%	79–100%		NPV	97%	89–100%	
NPV	96%	87–99%		NPV							

Explanation of symbols • Symbolforklaring • Interpretation der Symbole • Significado de los símbolos • Symbolien selitykset • Explication des symboles •
Symbolforklaring • Wyjaśnienie symboli • Teckenförklaringar

Use By
Anvendes inden
Verwendbar bis
Fecha de caducidad
Viimeinen käyttöpäivä
Utiliser jusque
Brukes før
Użyć przed
Använd före

Temperature limitation
Temperaturbegränsning
Temperaturbereich
Límite de temperatura
Lämpötilarajoitus
Limites de température
Temperaturbegrensning
Przestrzegać zakresu temperatury
Temperaturbegränsning

Manufacturer
Producent
Producent
Fabricante
Valmistaja
Fabricant
Produsent
Produsent
Tillverkare

Batch code
Batchnummer
Chargenbezeichnung
Codigo de lote
Erå^a
Code du lot
Parti
Kod partii
Lot nummer

Catalogue number
Katalognummer
Bestellnummer
Número de catálogo
Tuotenumero
Référence du catalogue
Katalognummer
Numer katalogowy
Katalognummer

Swab
Podepind
Tupfer
Hisopo
Näytteenottotikku
Écouvillon
Vattpinne
Wymaz
Provtagningspinne

Extraction buffer
Extraktionspuffer
Tampón de extracción
Uuttopuskuri
Tampon d'extraction
Tampone di estrazione
Bufor ekstrakcyjny
Tampão de Extracção
Extraktionsbuffert

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© Actim Oy. All rights reserved.
Actim is a trademark of Actim Oy.

AOACE36031-1 2021-10-28

CE