
INSTRUCTIONS FOR USE


Actim® Partus

EN Instructions for use	3	LV Lietošanas pamācība	68
CS Návod k použití	8	NL Gebruiksaanwijzing	73
DA Brugervejledning	13	NO Bruksanvisning	79
DE Gebrauchsanweisung	19	PL Instrukcja użytkowania	84
EL Οδηγίες χρήσης	24	PT Instruções de utilização	90
ES Instrucciones de uso	30	RO Instructiuni de utilizare	95
ET Kasutusjuhend	36	RU Инструкция по применению	101
FI Käyttöohje	41	SK Návod na použitie	107
FR Mode d'emploi	46	SL Navodila za uporabo	112
HU Használati utasítás	51	SV Bruksanvisning	117
IT Istruzioni per l'uso	57	TR Kullanım talimatları	122
LT Naudojimo instrukcija	62		

STRUCTURE OF THE DIPSTICK

TEST PROCEDURE AND RESULTS

INSTRUCTIONS FOR USE

The numbers and letters refer to the pictures on the inside cover.

STRUCTURE OF THE DIPSTICK

[A] Dip area [B] Result area [C] Test line [D] Control line [E] Handling area

INTENDED USE

Actim® Partus is a visually interpreted, manual qualitative immunochromatographic rapid test for the detection of phosphorylated IGFBP-1 (Insulin-like growth factor-binding protein 1) in cervical secretions during pregnancy when the fetal membranes are intact. The test is intended for professional use in clinical laboratories and near-patient testing (NPT) settings to help assess the risk of preterm or imminent birth starting from 22 weeks of gestation when signs and symptoms of preterm birth occur. **For *in vitro* diagnostic use only.**

BACKGROUND

Preterm birth (birth before the 37th week of gestation) is responsible for most of the non-chromosomal perinatal morbidity and mortality¹. The prediction of preterm birth remains a great challenge in obstetrics. Decidualized endometrial stromal cells synthesize the phosphorylated forms of IGFBP-1 (phIGFBP-1) while amniotic fluid contains substantial quantities of non- and less phosphorylated forms of IGFBP-1²⁻⁶. When delivery is approaching, the fetal membranes begin to detach from the decidua parietalis, and small amounts of phIGFBP-1 begin to leak into cervical secretions. In Actim Partus test a cervical specimen is taken with a sterile polyester swab during sterile speculum examination and extracted into the extraction buffer. The presence of phIGFBP-1 in the solution is detected using a dipstick.

PRINCIPLE OF THE TEST

The test is based on immunochromatography. It involves two monoclonal antibodies to human phIGFBP-1. One is bound to blue latex particles (the detecting antibody). The other is immobilized on a carrier membrane (the capture antibody) to catch the complex of antigen and latex-labeled antibody and indicate a positive

result. When the dip area of the dipstick is placed in an extracted specimen, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains phIGFBP-1, it binds to the antibody bound to latex particles. The particles are carried by the liquid flow and, if phIGFBP-1 is bound to them, they bind to the capture antibody. A blue line (test line) will appear in the result area, if the concentration of phIGFBP-1 in the sample exceeds the cut-off of the test. A blue control line confirms successful completion of the test.

KIT COMPONENTS

Actim Partus kit contains 10 test packs (Cat. no 31931ETAC) with instructions for use.

The components for performing one test are packed in one Actim Partus test pack (31921ETAC):

- **SWAB** 1 sterile polyester swab for specimen collection (EO).
- **EXTR BUF** 1 tube of extraction buffer (0.5 ml). This phosphate-buffered solution contains bovine serum albumin (BSA), protease inhibitors and preservatives (please see Warnings and precautions section).
- **DIP** 1 dipstick in a sealed aluminum foil pouch with desiccant.

Materials needed but not provided: Timer

STORAGE

Store the test kit at +2...+25 °C. When stored unopened, each component can be used until the expiry date marked on both the test pack and the kit box. The test kit and the test packs can also be stored for 2 months at +25...+30 °C, as long as the expiration date is not exceeded. All reagents are ready to use. Let the reagents reach room temperature before use, if stored at +2...+8 °C. Use the dipsticks shortly after their removal from the aluminum foil pouch. Do not freeze.

Before performing a test, read and follow the test instructions for use.

SPECIMEN COLLECTION AND EXTRACTION

If the kit components are stored refrigerated, allow the aluminum foil pouch and the extraction buffer to reach room temperature (+15...+25 °C) before use. The specimen is cervical secretion that is extracted into the provided extraction buffer. The specimen should be collected prior to performing digital examination and/

or transvaginal ultrasound. Take care not to touch anything with the swab before taking the sample.

1. Take the cervical secretion specimen using a sterile polyester swab (provided in the kit) from the cervical os during a sterile speculum examination. The swab should be left in the **external os for 10-15 seconds** to allow it to absorb the secretion specimen.
2. Open the extraction buffer tube and extract the specimen immediately from the swab by swirling the swab vigorously **in the extraction buffer for 10-15 seconds**. Press the swab against the wall of the extraction buffer tube to remove any remaining liquid from the swab. Discard the swab after extraction.
NOTE! Do not leave the swab to the tube.

SPECIMEN STORAGE

Specimens should be tested as soon as possible after they have been extracted, and no later than 4 hours following the specimen extraction. They should be stored at +2...+25 °C. If a specimen cannot be tested within this time it should be frozen. The frozen specimens are stable for 24 months (at or below -20°C). After thawing, the specimens should be mixed and tested as described below.

TEST PROCEDURE AND INTERPRETATION OF THE RESULTS

3. Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper purple part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch. Place the yellow dip area into the extracted specimen.
4. Hold the dipstick **there until you see the liquid front enter the result area**. Remove the dipstick from the solution and place it horizontally on a clean surface.
5. The result can be interpreted as positive as soon as two blue lines become visible in the result area. Negative result should be read at **5 minutes**. **Do not pay attention to any lines appearing later than 5 minutes**. If **two blue lines, the test line and control line**, appear, the test result is **positive (5A)**. Any faint but continuous line should be interpreted as positive. If **one blue line, the control line**, appears, the test result is **negative (5A)**. If **the control line does not appear**, the test is **invalid (5B)**.

QUALITY CONTROL

Regular use of Actim® Partus Controls (Cat. No. 31900ETAC) is recommended.

Actim Partus Controls kit includes 3 freeze-dried controls (a negative, a low positive, and a high positive control) and one reconstitution solution tube. The controls are provided separately.

Please consult Actim Partus Controls instructions for use for further information and instructions on how to use the controls.

Quality control shall be performed according to the applicable national or local regulations or accreditation requirements and laboratory standard Quality Controls Procedures. It is recommended to perform a control test at least with new reagent lot and shipment.

LIMITATIONS OF THE TEST

- The test is intended for *in vitro* diagnostic use only.
- For professional use only.
- No quantitative interpretation should be made based on the test results.
- Before performing the test ensure that the fetal membranes are intact (for example with Actim PROM test), because with ruptured fetal membranes Actim Partus test will also give a positive result.
- Patients with moderate or heavy vaginal bleeding should not be tested. It is recommended that a sample be taken when bleeding has stopped, and the extract is essentially blood-free.
- The test result indicates the risk at the time of specimen collection and changes in the patient's condition may later affect the final outcome of the pregnancy.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.

WARNINGS AND PRECAUTIONS

- Do not eat, drink or smoke in the area where the kit and specimens are handled.
- Wear protective personal equipment when handling and performing a test. Wash hands after specimen and test handling.
- The extraction buffer contains <0.01 % reaction mass of: 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) [EC No: 611-341-5]

(Skin Sens 1, Aquatic Chronic 3). May cause an allergic skin reaction (H317). Harmful to aquatic life (H402). Harmful to aquatic life with long lasting effects (H412). Avoid breathing mist/vapour (P261). Contaminated work clothing should not be allowed out of the workplace (P272). Wear eye protection/protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents/container in accordance with local regulation (P501).

Warning

- Do not use the product after the expiry date marked on the test pack or kit box. Do not exceed the stated stability periods for opened components.
- Do not mix the components from different tests and lot numbers. Components are disposable, do not reuse them.
- The test requires about 150 µl of extracted sample to ensure proper performance of the test.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing.
- Do not use the dipstick if its aluminum foil pouch or the seals of the pouch are not intact. Do not use the swab if the package is damaged.
- Use only the swab provided with the kit.
- Improper sampling and/or testing procedure may lead to false results.
- Care must be taken when placing the dipstick in the extraction buffer tube. The upper part of the dipstick must stay dry.
- If the control line does not appear, the test is invalid, and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly it is recommended that the test be repeated.
- Do not pay attention to any lines appearing after 5 minutes.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

- Temperatures above or below the recommendations can cause inaccurate results.
- In case of serious incident, please report it to the manufacturer or its representative and/or national authority.

PERFORMANCE OF THE TEST

Analytical Sensitivity

The analytical sensitivity of the Actim Partus test was identified by evaluating different concentrations of phIGFBP-1 in extracted samples on three different lots of the Actim Partus test. Two different operators each interpreted ten devices run at each concentration under various lighting conditions for a total of 60 determinations per level. The cut-off of Actim Partus test is approximately 10 µg/l in the extracted specimen.

Measuring range

The measuring range of the Actim Partus is approximately 10-10 000 µg/l in the extracted specimen.

Analytical Specificity

Analytical specificity (cross-reactivity) was tested with human IGFBP proteins at a concentration level of 5000 µg/l of IGFBP-2, -4, -5 and -6 proteins, and at 50 000 µg/l of IGFBP-3 protein using one lot of the Actim Partus test. No cross-reactivity was seen using human IGFBP-2, -3, -4, -5 and -6 proteins. The Actim Partus test is specific to human IGFBP-1.

Repeatability and reproducibility

Repeatability (intra-assay) was evaluated with a panel of specimens at different concentration levels of phIGFBP-1. Specimens were tested as 10 replicates during the same day, using three different lots of Actim Partus. Repeatability of negative and positive specimens was 100%, and repeatability of specimens at borderline concentration was 57%.

Reproducibility (inter-assay) study of the Actim Partus test was conducted at three different sites. Test operators (n=9) tested the same panel of samples on five different days using one lot of Actim Partus test. A total of 630 tests were performed (210 per site) with a total of 90 tests per specimen. The overall

reproducibility of Actim Partus test is 97.5% (614/630) with no significant differences within runs (replicates tested by one operator), between runs (5 different days), between sites (3 sites) or between operators (9 operators).

Diagnostic performance

Actim Partus test has been evaluated in several clinical studies. Selections of typical studies are shown in table 1. In these studies, Actim Partus test was used in predicting preterm birth within 7 days and in predicting labor before 32–37 weeks (table 1.).

Interference testing

Take care not to touch anything with the swab before taking the specimen. Do not contaminate the swab or cervicovaginal specimen with lubricants or creams, as they may physically interfere with absorption of the specimen onto the swab and/or affect the test performance.

TABLE 1.

Study	GA (Gestational age)	Outcome	Sensitivity (%)	Specificity (%)	Positive predictive value (PPV) (%)	Negative predictive value (NPV) (%)
Azlin et al. 2010 ⁷	24-36	7 days	80	94	57	98
Brik Spinelli et al. 2010 ⁸	24-34	7 days	73	66	22	95
		< 32 weeks	76	66	18	96
Tanir et al. 2009 ⁹	24-37	7 days	93	79	56	98
		< 34 weeks	70	75	48	89
Eroglu et al. 2007 ¹⁰	24-35	7 days	83	84	42	97
		< 35 weeks	70	88	58	92
Akercan et al. 2004 ¹¹	24-36	< 37 weeks	78	87	73	90
Altinkaya et al. 2009 ¹²	24-34	< 37 weeks	70	87	56	93
Ting et al. 2007 ¹³	24-34	7 days	69	78	39	92
Lembet et al. 2002 ¹⁴	20-36	7 days	94	85	83	94
		< 37 weeks	90	94	94	89

The following substances, conditions and micro-organisms were tested with the Actim Partus test and were found not to affect the Actim Partus test performance, when tested at the concentrations shown in Table 2.

TRACEABILITY

A primary phIGFBP-1 antigen preparate was purified from the human decidua tissue and the concentration was determined by the UV spectrophotometric analysis using wavelength 280 nm. This preparate was used to prepare the primary reference standards for the quantitative phIGFBP-1 ELISA method to which the Actim Partus is traceable to.

DISPOSAL

Dispose of the contents according the national and local law.

All patient specimens and used swabs, dipsticks and extraction buffer tubes should be handled and disposed of as potentially infectious material.

TABLE 2.

Substance/Condition/ Micro-organism**Tested concentration**

Semen	100 % v/v
Urine	100 % v/v
Surgical lubricant water-based (Aquasoniq 100, Optilube, Surgilube, HR Lubricating Jelly)	6 % w/v
Betadine	50 % w/v
Skin oil (Ceridal Lipolotion)	50 % v/v
Personal lubricant water-based (K-Y Jelly, Klick)	25 % w/v
Personal lubricant silicone-based (Magic Glide)	6 % w/v
Hexochlorophene	250 mg/ml
Hexochlorophene cream (Hibitane)	100 mg/ml
Econazole cream (Gyno-Pevaryl)	250 mg/ml
Prostaglandin E2	100 g/ml
Ampicillin	100 µg/ml
Erythromycin	60 µg/ml
Dexamethasone	200 µg/ml
Betamethasone	575 ng/ml
Naproxen	500 µg/ml
Nifedipine	5 µg/ml
Magnesium sulfate	50 µg/ml
Oxytocin	100 Units/ml
Terbutaline	1 mg/ml
17-alpha hydroxyprogesterone caproate	175 ng/ml
Progesterone capsule oral use	75 ng/ml
Progesterone capsule vaginal use	250mg/ml
Indomethacin	36 µg/ml
Bilirubin	500 µg/ml
Shower and bath products (Rexona, Dove, Palmolive, Adidas, Daily)	1 % w/v

Substance/Condition/ Micro-organism**Tested concentration**

Baby Oil (Natusan)	50 % v/v
Baby Powder (Natusan)	50 % w/v
Feminine deodorant (Vagi-Gard)	50% w/v
Feminine deodorant film (VCF Dissolving Feminine Deodorant Film)	25 % w/v
Vaginal film (VCF Dissolving Vaginal Cleansing Film)	50% w/v
Vaginal gel (Refresh)	2 % w/v
Candida albicans	11.2* 10 ⁸ CFU/ml
Gardnerella vaginalis	11.8 * 10 ⁸ CFU/ml
Neisseria gonorrhoea	8.6 *10 ⁸ CFU/ml
Chlamydia trachomatis	5.0 * 10 ³ TCID ⁵⁰
HSV-1	5.0 * 10 ³ TCID ⁵⁰
HSV-2	5.0 * 10 ³ TCID ⁵⁰

Materials used in Actim Partus kit and its components:

Paper: instructions for use.

Cardboard: Kit box.

Plastic: test pack pouch, extraction buffer tube.

Several (not to be recycled): unused dipstick, unopened extraction buffer tube, swabs, dipstick pouch.

When used in accordance with Good Laboratory Practice, good occupational hygiene and instructions for use, the reagents should not present a hazard to health.

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TESTOVANÁ KONCENTRACE

Indomethacin	36 µg /ml
Bilirubin	500 µg /ml
Shower and bath products (Rexona, Dove, Palmolive, Adidas, Daily)	1 % w/v
Baby Oil (Natusan)	50 % v/v
Baby Powder (Natusan)	50 % w/v
Feminine deodorant (Vagi-Gard)	50% w/v
Feminine deodorant film (VCF Dissolving Feminine Deodorant Film)	25 % w/v
Vaginal film (VCF Dissolving Vaginal Cleansing Film)	50% w/v
Vaginal gel (Refresh)	2 % w/v
Candida albicans	11.2* 10 ⁸ CFU/ml
Gardnerella vaginalis	11.8 * 10 ⁸ CFU/ml
Neisseria gonorrhoea	8.6 *10 ⁸ CFU/ml
Chlamydia trachomatis	5.0 * 10 ³ TCID ⁵⁰
HSV-1	5.0 * 10 ³ TCID ⁵⁰
HSV-2	5.0 * 10 ³ TCID ⁵⁰

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1. Born too soon: decade of action on preterm birth. Geneva: World Health Organization; 2023. Licence: CC BY-NC-SA 3.0 IGO. (<https://creativecommons.org/licenses/by-nc-sa/3.0/igo/>)
2. Jones et al., 1991. Jones JI, D'Ercole AJ, Camacho-Hubner C and Clemons DR (1991) Phosphorylation of insulin-like growth factor (IGF) binding protein 1 in cell culture and *in vivo*: effects on affinity for IGF-I Proceedings of the National Academy of Sciences USA 88 7481–7485
3. Westwood et al., 1994. Westwood M, Gibson JM, Davies AJ, Young RJ and White A (1994) The phosphorylation pattern of insulin-like growth factor-binding protein 1 in normal plasma is different from that in amniotic fluid and changes during pregnancy Journal of Clinical Endocrinology and Metabolism 79 1735–1741
4. Frost and Tseng, 1991. Frost RA and Tseng L (1991) Insulin-like growth factor-binding protein 1 is phosphorylated by cultured human endometrial stromal cells and multiple protein kinases *in vitro*. Journal of Biological Chemistry 266 18 082–18 088
5. Martina et al., 1997 Martina NA, Kim E, Chitkara U, Wathen NC, Chard T and Giudice LC (1997) Gestational age-dependent expression of insulin-like growth factorbinding protein 1 (PHIGFBP-1) phosphoisoforms in human extraembryonic cavities, maternal serum, and decidua suggests decidua as the primary source of PHIGFBP-1 in these fluids during early pregnancy Journal of Clinical Endocrinology and Metabolism 82 1894–1898
6. Westwood et al., 1998 Westwood M, Gibson JM, Aplin JD and White A (1998) PHIGFBP-1 regulation of IGF-I bioavailability in pregnancy Journal of Reproduction and Fertility Abstract Series 21 Abstract 1
7. Nor Azlin MI et al. Role of pHIIGFBP-1 and ultrasound cervical length in predicting pre-term labour. Journal of Obstetrics and Gynaecology (2010) 30:456-460.
8. Brik Spinelli M et al. Phosphorylated insulin-like growth factor binding protein-1 and cervical measurement in women with threatening preterm birth. Acta Obstet Gynecol Scand (2010) 89:268-74.
9. Tanir HM, Sener T, Yildiz Z. Cervical phosphorylated insulin-like growth factor binding protein-1 for the prediction of preterm delivery in symptomatic cases with intact membranes. J Obstet Gynaecol Res (2009) 1:66-72.
10. Eroglu D et al. Prediction of preterm delivery among women with threatened preterm labor. Gynecol Obstet Invest (2007) 64:109-116.
11. Akercan F et al. Value of cervical phosphorylated insulinlike growth factor binding protein-1 in the prediction of preterm labor. J Reprod Med (2004) 49:368-372.
12. Altinkaya O et al. Cervical phosphorylated insulin-like growth factor binding protein-1 in prediction of preterm delivery. Arch Gynecol Obstet (2009) 279:279-283. Epub 2008 Jun 13.
13. Ting HS et al. Comparison of bedside test kits for prediction of preterm delivery: phosphorylated insulin-like growth factor binding protein-1 (pIGFBP-1) test and fetal fibronectin test. Ann Acad Med Singapore (2007) 36:399-402.
14. Lembet A et al. New rapid bed-side test to predict preterm delivery: phosphorylated insulin-like growth factor binding protein-1 in cervical secretions. Acta Obstet Gynecol Scand (2002) 81:706-712.

Explanation of symbols • Sümbole selgitused • Vysvětlení symbolů • Symbolforklaring • Erklärung der Symbole • Ενεργήση των συμβόλων • Explicación de los símbolos • Symbolien selitykset • Explication des symboles • Szimbólumok magyarázata • Spiegazione dei simboli • Simboli aiškinimas • Simbolu skaidrojums • Uitleg van symbolen • Symbolforklaringer • Objasnenie symboli • Explicação dos símbolos • Explicația simbolurilor • Описание символов • Vysvetlivky symbolov • Razlagi simbolov • Symbolförförklaring • Sembollerin açıklamaları

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Contains <0.01 % reaction mass of: 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) [EC No: 611-341-5]

Sisaldb < 0,01 massiprotsenti järgmist reaktivi: 5-kloro-2-metüül-2H-isotiasool-3-oon ja 2-metüül-2H-isotiasool-3-oon (3 : 1) [EÜ nr: 611-341-5]

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© Actim Oy. All rights reserved.
Actim® is a trademark of Actim Oy.
Actim Oy