

Actim® CRP

STRUCTURE OF DIPSTICK

Test lines [3] [4] Control line

TEST PROCEDURE AND RESULTS

INSTRUCTIONS FOR USE

Numbers [1]- [10] refer to illustrations.

STRUCTURE OF DIPSTICK

[1] Dip area [2] Result area [3] Test lines [4] Control line

The dipstick is used as a visually interpretive, semiquantitative immunochromatographic dipstick test, which is used for detection and monitoring of CRP concentrations in whole blood samples. The test is intended for professional use.

KIT COMPONENTS

The Actim CRP kit (31031ETAC) contains 20 CRP test packs with instructions for use. The components of each Actim CRP test pack (31021ETAC) are:

- 1x Actim CRP dipstick in sealed aluminium pouch with desiccant.
- 1x [DIL] [BUF] One tube of dilution buffer (0.5 mL). The buffer contains buffer salts, bovine serum albumin (BSA) and preservatives (please see Warnings and precautions section).
- 1x [CAPL] Each kit contains a vial of 10 µl end-to-end capillaries (25). Alternatively, other 10 µl capillaries or a pipette may be used for sample collection (not available with the kit).

STORAGE

Store the test kit at +2...+25 °C. Stored unopened, each component can be used until the expiry date marked on the container. After the pouch is opened, the dipstick from the pouch just prior to use.

SPECIMEN COLLECTION

The sample is 10 µl of whole blood diluted in dilution buffer. The diluted sample can be kept at room temperature for a hours.

Sampling from a fingertip with the end-to-end capillary
Use a lancet or a bandage device to obtain a drop of blood from a fingertip. Take the sample with the 10 µl capillary [5]. Make sure that the capillary becomes filled with blood. Immediately drop the capillary into the dilution buffer tube. Mix the sample into the buffer by inverting the tube upside down 10-15 times until the capillary is totally empty. When the CRP concentration is above 80 mg/L, the sample is above the capillary [6].

Sampling from a whole blood sample tube
EDTA, citrate or heparin blood can also be used as a sample. The sample can be pipetted from the blood sample into a 10 µl pipette. Dispense the sample into the dilution buffer tube and mix the sample into the buffer.

- If stored refrigerated, allow the aluminium foil pouch and the buffer tube to reach room temperature. Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper green part of the dipstick. The dipstick must be used shortly after opening the pouch.
- Allow the specimen to reach room temperature before testing. Place the yellow dip area (the lower end of the dipstick) into the diluted sample [7] and hold it there **until you see the liquid front enter the result area [8]**. Remove the dipstick from the solution and place it in a horizontal position [8]. The result is available in **5 minutes [9]**. **Do not pay attention to any lines appearing later than 5 minutes [10]**.

If, in addition to the red control line there are

- no blue lines**, the serum CRP concentration of the sample is **<10 mg/L**
- one blue line**, the serum CRP concentration is **10-40 mg/L**
- two blue lines**, the serum CRP concentration is **40-80 mg/L**
- three blue lines**, the serum CRP concentration is **>80 mg/L**

In the control line *does not appear*, the test is invalid [10]

The intensities of the blue lines vary. When the CRP concentration is low (10-15 mg/L) only one first blue line appears in 5 minutes. As the concentration of CRP increases, the intensity of the first line increases. A second blue line begins to appear when the concentration exceeds 40 mg/L. Corresponding concentrations between 40 and 80 mg/L are not clearly distinguishable. A third blue line begins to appear when the CRP concentration exceeds 80 mg/L.

LIMITATIONS OF THE TEST

- The test is intended for *in vitro* diagnostic use only.
- For professional use only.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.

WARNINGS AND PRECAUTIONS

- Do not eat, drink or smoke in the area where the kit and specimens are handled.
- Wear protective personal equipment when handling and performing a test. Wash hands after specimen and test handling.
- The dilution buffer contains <0.01% reaction mass of 5-β-Chloro-2-methyl-2H-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one [31] [EC No. 611-341-5] (Skin Sens 1) (Aquatic Chronic 3). May cause an allergic skin reaction (H317). Harmful to aquatic life (H412). Harmful to aqmatic life with long lasting effects (H412). Avoid breathing mist/vapour (P273). Wear eye protection/protective clothing should be allowed out of the workplace (P272). Wear eye protection/protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents/container in accordance with local regulation (P501).
- Do not use the product after the expiry date marked on the test pack or kit box. Do not exceed the stated stability periods for opened components.
- Do not mix the components from different tests and lot numbers. Components are disposable, do not reuse them.
- Carry enough 150 µl of diluted sample to ensure proper performance of the test.
- Do not care be taken when placing the dipstick in the sample tube. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice any condensation on the dipstick area before testing.
- Do not use a dipstick if its aluminium foil pouch or the seals are not intact.
- The red control lines is in the upper part of the result area. Depending on the CRP concentration, up to three blue lines indicating different concentrations of CRP become visible starting from the lower part of the result area. Appearance of a control line confirms the result.
- If a control line does not appear the test is invalid, and should be repeated using another dipstick.
- The intensities of the test lines may vary considerably. Please note that the second and the third line never become as intense as the first line.
- Sometimes erythrocytes may cause red coloring in the lower part of the result area. This is normal and does not have an effect on the performance of the test.
- If the test result cannot be interpreted clearly (e.g. if the lines are blotched or uneven) it is recommended that the test be repeated.
- The result of a blood sample is available 5 minutes after 5 minutes have elapsed. However, do not pay attention to any lines appearing after 5 minutes.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

PRINCIPLE OF THE TEST
CRP is an acute phase protein produced rapidly in response to inflammation, infection and tissue damage. It is used in particular to differentiate bacterial infections from virus infections. The test is based on immunochromatography. It involves monoclonal antibodies to human CRP. One is bound to blue latex particles (the detecting label). In the carrier membrane there are three CRP-specific antibody zones, to which the blue latex particles will bind. The particles of CRP, when the top area of the dipstick is placed in a diluted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains CRP, it binds to the antibody labeled with latex particles. The particles are carried by the liquid flow and, if CRP is bound to them, they bind to the catching antibody zones. The more CRP the sample contains, the more blue lines become visible.

The red control line confirms the correct performance of the test.

PERFORMANCE OF THE TEST

Analytical sensitivity and measuring range
The lowest detectable amount of CRP in serum samples is 10 mg/L. One blue line appears when the concentration is 10 mg/L. The second blue line appears when the concentration is 40 mg/L. The third line becomes visible if the serum CRP concentration is 80 mg/L. Hook effect has been tested with CRP-spiked samples. Samples containing up to 750 mg/L give consistently the result >80 mg/L. CRP test is traceable to the reference material CRM-DA470.

Repeatability and reproducibility
Intra-laboratory reproducibility
Inter- (reproducibility) assay precision studies have confirmed that the three different cut-off values of the CRP test are approximately 10, 40 and 80 mg/L.

Interference testing
No interferences of bilirubin, rheumatoid factor (RF), lipid and SAP (Serum Amyloid P) have been observed.

Diagnostic performance
Actim CRP test was evaluated at three different hospitals in Finland. At hospitals 1 and 2 EDTA blood samples were tested against a quantitative method. The evaluations showed that the Actim CRP results are comparable to the quantitative measurement. Results are shown in FIG 1 and 2. At hospital 3 Actim CRP test was performed from fingertip blood and EDTA blood samples. The evaluation showed that the Actim CRP results are comparable to the quantitative measurement. Results are shown in FIG 3 and 4.

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- Wear protective personal equipment when handling and performing a test. Wash hands after specimen and test handling.
- The dilution buffer contains <0.01% reaction mass of 5-β-Chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one [31] [EC No. 611-341-5] (Skin Sens 1) (Aquatic Chronic 3). Can cause an allergic skin reaction (H317). Harmful to aquatic life (H412). Harmful to aquatic life with long lasting effects (H412). Avoid breathing mist/vapour (P273). Wear eye protection/protective clothing should be allowed out of the workplace (P272). Wear eye protection/protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents/container in accordance with local regulation (P501).
- Do not use the product after the expiry date marked on the test pack or kit box. Do not exceed the stated stability periods for opened components.
- Do not mix the components from different tests and lot numbers. Components are disposable, do not reuse them.
- Carry enough 150 µl of diluted sample to ensure proper performance of the test.
- Do not care be taken when placing the dipstick in the sample tube. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice any condensation on the dipstick area before testing.
- Do not use a dipstick if its aluminium foil pouch or the seals are not intact.
- The red control lines is in the upper part of the result area. Depending on the CRP concentration, up to three blue lines indicating different concentrations of CRP become visible starting from the lower part of the result area. Appearance of a control line confirms the result.
- If a control line does not appear the test is invalid, and should be repeated using another dipstick.
- The intensities of the test lines may vary considerably. Please note that the second and the third line never become as intense as the first line.
- Sometimes erythrocytes may cause red coloring in the lower part of the result area. This is normal and does not have an effect on the performance of the test.
- If the test result cannot be interpreted clearly (e.g. if the lines are blotched or uneven) it is recommended that the test be repeated.
- The result of a blood sample is available 5 minutes after 5 minutes have elapsed. However, do not pay attention to any lines appearing after 5 minutes.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

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