
INSTRUCTIONS FOR USE


Actim[®]PROM

| | | | | | |
|----|-----------------------|----|----|---------------------------|-----|
| EN | Instructions for use | 3 | LT | Naudojimo instrukcijos | 60 |
| ET | Kasutusjuhend | 7 | LV | Lietošanas pamācība | 65 |
| CS | Návod k použití | 12 | NL | Gebruiksaanwijzing | 70 |
| DA | Brugervejledning | 17 | NO | Bruksanvisning | 75 |
| DE | Gebrauchsanweisung | 21 | PL | Instrukcja użytkowania | 79 |
| EL | Οδηγίες χρήσης | 26 | PT | Instruções de utilização | 84 |
| ES | Instrucciones de uso | 31 | RO | Instrucțiuni de utilizare | 89 |
| FI | Käyttöohje | 36 | RU | Инструкция по применению | 94 |
| FR | Mode d'emploi | 41 | SK | Návod na použitie | 100 |
| HU | Használati utasítás | 46 | SL | Navodila za uporabo | 105 |
| IT | Istruzioni per l'uso | 51 | SV | Bruksanvisning | 109 |
| IS | Notkunarleiðbeiningar | 56 | TR | Kullanım talimatları | 114 |

STRUCTURE OF THE DIPSTICK

TEST PROCEDURE AND RESULTS

INSTRUCTIONS FOR USE

The numbers and letters refer to the pictures on the inside cover.

STRUCTURE OF THE DIPSTICK

[A] Dip area [B] Result area [C] Test line [D] Control line [E] Handling area

INTENDED USE

Actim® PROM is a visually interpreted, manual qualitative immunochromatographic rapid test for detection of amniotic fluid in cervicovaginal secretions during pregnancy. Actim PROM detects IGFBP-1 (Insulin-like growth factor-binding protein 1), which is a major protein in amniotic fluid and a marker of the amniotic fluid in cervicovaginal secretions. The test is intended for professional use in clinical laboratories and near-patient testing (NPT) settings to help diagnose the rupture of fetal membranes (ROM) during pregnancy when a rupture of fetal membranes is suspected. **For *in vitro* diagnostic use only.**

BACKGROUND

Premature rupture of fetal membranes (PROM) is defined as a rupture of fetal membranes before onset of the labor. PROM is associated with increased risk of intra-uterine infection and thus it increases the risk of both maternal and perinatal morbidity and mortality. If the rupture of membranes occurs before 37 weeks of gestation, it is referred to as preterm PROM (PPROM)¹. PROM occurs in approximately 5-10% of all pregnancies and PPROM in approximately 1-3%. Management of a patient with PROM or PPROM requires evaluation of the relative risks of preterm birth and prolongation of gestation. The Actim PROM test is intended to assist the healthcare professional in making an informed decision about how to manage the patient.

The concentration of IGFBP-1 in amniotic fluid is 100–1,000 times higher than in maternal serum². Seminal fluid and urine do not contain a measurable amount of IGFBP-1, and it is not present in the cervicovaginal secretions if the fetal membranes are intact². After the rupture of fetal membranes, amniotic fluid with a high concentration of IGFBP-1 mixes with cervicovaginal fluids². In the Actim PROM test, a specimen is taken with a sterile polyester swab from cervicovaginal fluid, and the specimen is extracted into the extraction buffer. The presence of IGFBP-1 in the extracted specimen is detected using a dipstick.

PRINCIPLE OF THE TEST

The test is based on immunochromatography. It involves two monoclonal antibodies to human IGFBP-1. One is bound to blue latex particles (the detecting antibody). The other is immobilized on a carrier membrane (the capture antibody) to catch the complex of antigen and latex-labeled antibody and indicate a positive result. When the dip area of the dipstick is placed in an extracted specimen, the dipstick absorbs liquid, which starts to flow up the dipstick. If the specimen contains IGFBP-1, it binds to the antibody bound to latex particles. The particles are carried by the liquid flow and, if IGFBP-1 is bound to them, they bind to the capture antibody. A blue line (test line) will appear in the result area, if the concentration of IGFBP-1 in the specimen exceeds the cut-off of the test. A blue control line confirms successful completion of the test.

KIT COMPONENTS

Actim PROM kit contains 20 (Cat. No 30832ETAC) or 10 test packs (Cat. No. 30831ETAC) with instructions for use.

The components for performing one test are packed in one Actim PROM test pack (Cat. No 30821ETAC):

- **[SWAB]** 1 sterile polyester swab for specimen collection (EO).
- **[EXTR]** **[BUF]** 1 tube of extraction buffer (0.5 ml). This phosphate-buffered solution contains bovine serum albumin (BSA), protease inhibitors and preservatives (please see Warnings and precautions section).
- **[DIP]** 1 dipstick in a sealed aluminum foil pouch with desiccant.

Materials needed but not provided: Timer

STORAGE

Store the test kit at +2...+25 °C. When stored unopened, each component can be used until the expiry date marked on both the test pack and the kit box. The test kit and the test packs can also be stored for 2 months at +25...+30 °C, as long as the expiration date is not exceeded. All reagents are ready to use. Let the reagents reach room temperature before use, if stored at +2...+8 °C. Use the dipsticks shortly after their removal from the aluminum foil pouch. Do not freeze.

Before performing a test, read and follow the test instructions for use.

SPECIMEN COLLECTION AND EXTRACTION

If the kit components are stored refrigerated, allow the aluminum foil pouch and the extraction buffer to reach room temperature (+15...+25 °C) before use. The specimen is a cervicovaginal secretion that is extracted into the provided extraction buffer. A specimen is collected using a sterile polyester swab (provided in the kit). The specimen should be collected prior to performing digital examination and/or transvaginal ultrasound. Take care not to touch anything with the swab before taking the specimen.

1. Separate the labia and carefully insert the tip of the swab into the vagina toward the posterior fornix until resistance is met. Alternatively, the specimen can be taken from the posterior fornix during a sterile speculum examination. The swab should be held in **the posterior fornix for 10–15 seconds** to allow it to absorb the cervicovaginal secretion.
2. Open the extraction buffer tube and extract the specimen immediately from the swab by swirling the swab vigorously in **the extraction buffer for 10–15 seconds**. Press the swab against the wall of the extraction buffer tube to remove any remaining liquid from the swab. Discard the swab after extraction. **NOTE!** Do not leave the swab in the tube.

SPECIMEN STORAGE

Specimens should be tested as soon as possible after they have been extracted, and no later than 4 hours following the specimen extraction. They should be stored at +2...+25 °C. If a specimen cannot be tested within this time, it should be frozen. The frozen specimens are stable for 24 months (at or below -20°C). After thawing, the specimens should be mixed and tested as described below.

TEST PROCEDURE AND INTERPRETATION OF THE RESULTS

3. Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper turquoise part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch. Place the yellow dip area into the extracted specimen.
4. Hold the dipstick in the tube **until you see the liquid front enter the result area**. Remove the dipstick from the buffer and place it horizontally on a clean surface.

5. The result can be interpreted as positive as soon as two blue lines become visible in the result area. Negative result should be read at **5 minutes. Do not pay attention to any lines appearing later than 5 minutes. If two blue lines, the test line and the control line, appear, the test result is positive (5A).** Any faint but continuous line should be interpreted as positive. If **only the control line appears**, the test result is **negative (5A)** If **the control line does not appear**, the test is **invalid (5B)**.

QUALITY CONTROL

Regular use of Actim® PROM Controls (Cat. No. 30800ETAC) is recommended. Actim PROM Controls kit includes 3 freeze-dried controls (a negative, a low positive, and a high positive control) and one reconstitution solution. The controls are provided separately.

Please consult Actim PROM Controls instructions for use for further information and instructions on how to use the controls.

Quality control shall be performed according to the applicable national or local regulations or accreditation requirements and laboratory standard Quality Controls Procedures. It is recommended to perform a control test at least with new reagent lot and shipment.

LIMITATIONS OF THE TEST

- The test is intended for *in vitro* diagnostic use only.
- For professional use only.
- No quantitative interpretation should be made based on the test results.
- If rupture of fetal membranes has occurred **but the leakage of amniotic fluid has ceased** more than 12 hours before the specimen is taken, IGFBP-1 may have been degraded by proteases in the vagina and the test may give a negative result.
- A positive Actim PROM test result, although detecting the presence of amniotic fluid in the specimen, does not locate the site of the rupture.
- A negative test result is an indication on the present condition and cannot be used to predict the forthcoming.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.

WARNINGS AND PRECAUTIONS

- Do not eat, drink or smoke in the area where the kit and specimens are handled.
- Wear protective personal equipment when handling and performing a test. Wash hands after specimen and test handling.
- The extraction buffer contains <0.01% reaction mass of: 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) [EC No: 611-341-5] (Skin Sens 1, Aquatic Chronic 3). May cause an allergic skin reaction (H317). Harmful to aquatic life (H402). Harmful to aquatic life with long lasting effects (H412). Avoid breathing mist/vapour (P261). Contaminated work clothing should not be allowed out of the workplace (P272). Wear eye protection/protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents/container in accordance with local regulation (P501).
- Do not use the product after the expiry date marked on the test pack or kit box. Do not exceed the stated stability periods for opened components.
- Do not mix the components from different tests and lot numbers. Components are disposable, do not reuse them.
- The test requires about 150 µl of extracted specimen to ensure proper performance of the test.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick, if you notice a blue coloring in the result area before testing.
- Do not use the dipstick, if its aluminum foil pouch or the seals of the pouch are not intact. Do not use the swab if the package is damaged.
- Use only the swab that is provided with the kit.
- Improper sampling and/or testing procedure may lead to false results.
- Care must be taken when placing the dipstick in the extraction buffer tube. The upper part of the dipstick must stay dry.

- If the control line does not appear, the test is invalid, and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly, it is recommended that the test be repeated.
- Do not pay attention to any lines appearing after 5 minutes.
- All biological specimens and materials must be treated as potentially hazardous and disposed of in accordance with local authority guidelines.
- Temperatures above or below the recommendations can cause inaccurate results.
- In case of serious incident, please report it to the manufacturer or its representative and/or national authority.

PERFORMANCE OF THE TEST

Analytical Sensitivity

The analytical sensitivity of the Actim PROM test was identified by evaluating different concentrations of IGFBP-1 in extracted specimens on three different lots of the Actim PROM test. Two different operators each interpreted ten tests run at each concentration under various lighting conditions for a total of 60 determinations per level. The cut-off of Actim PROM test is approximately 25 µg/l in extracted specimen.

Measuring range

The measuring range of the Actim PROM is approximately 25-200 000 µg/l in extracted specimen.

Analytical Specificity

Analytical specificity (cross-reactivity) was tested with human IGFBP proteins at concentrations ranging from 10-5000 µg/l of each protein in extracted specimen using three different lots of the Actim PROM test. No cross-reactivity was seen using human IGFBP-2, -3, -4, -5 and -6 proteins. The Actim PROM test is specific to human IGFBP-1.

Repeatability and reproducibility

A panel of specimens consisting of samples of different IGFBP-1 concentration levels was evaluated for intra- (repeatability) and inter- (reproducibility) assay precision. For repeatability the samples were tested with 10 replicates during the same day using three different lots of the Actim PROM test. Repeatability of negative and positive specimens was 100% and repeatability of specimen at borderline concentration was 67 %. For reproducibility a study of the Actim PROM test was conducted at three separate sites. Test operators (n=9) tested each level on five different days using one lot of the Actim PROM test. A total of 360 tests were performed (120 per site) with a total of 45 tests per sample type. The overall reproducibility of the Actim PROM test is 97% (350/360) with no significant differences within runs (replicates tested by one operator), between runs (five different days), between sites (three sites) or between operators (nine operators).

Diagnostic performance

The Actim PROM test has been evaluated in several clinical studies. A selection of typical studies is shown in Table 1.

TABLE 1.

| | Sensitivity (%) | Specificity (%) | Positive predictive value (PPV) (%) | Negative predictive value (NPV) (%) |
|-----------------------------|------------------------|------------------------|--|--|
| Rutanen 1996 ³ | 100 | 95 | 93 | 100 |
| Ragosch 1996 ⁴ | 100 | 83 | 83 | 100 |
| Jain 1998 ⁵ | 100 | 89 | 76 | 100 |
| Kubota 1998 ⁶ | 95 | 93 | 90 | 96 |
| Ergemoglu 2004 ⁷ | 97 | 97 | 97 | 97 |

Interference testing

Take care not to touch anything with the swab before taking the specimen. Do not contaminate the swab or cervicovaginal specimen with lubricants or creams, as they may physically interfere with absorption of the specimen onto the swab and/or affect the test performance. Large amounts of meconium causing green

staining of the amniotic fluid may interfere with the detection of IGFBP-1 with Actim PROM.

The following substances, conditions and micro-organisms were tested with Actim PROM test and were found not to affect Actim PROM test performance, when tested at the concentrations shown in Table 2.

TABLE 2.

| Substance/Condition/ Micro-organism | Tested concentration |
|---|-------------------------------|
| Whole blood | 100% v/v |
| Semen | 100% v/v |
| Urine | 100% v/v |
| pH 3.5-8.5 | N/A |
| Surgical lubricant water-based (Aquasoniq 100, Optilube, Surgilube, HR Lubricating Jelly) | 50% w/v |
| Skin oil (Cerial Lipolotion) | 50% v/v |
| Pevaryl (active ingredient: econazol.nitras) | 30 mg/ml |
| Gyno-Trosyd (tioconazol) | 20 mg/ml |
| Flagyl (metronidazole) | 100 mg/ml |
| Canesten (clotrimazol) capsule | 40 mg/ml |
| Personal lubricant water-based (K-Y Jelly, Klick) | 25% w/v |
| Personal lubricant silicone-based (Magic Glide) | 3% w/v |
| Baby oil (Natusan) | 50% v/v |
| Baby powder (Natusan) | 50% w/v |
| Feminine deodorant (Vagi-Gard) | 50% w/v |
| Vaginal film (VCF Dissolving Vaginal Cleansing Film) | 0.1% w/v |
| Vaginal gel (RepHresh) | 25% w/v |
| Shower and bath products (Rexona, Dove, Palmolive, Adidas, Daily) | 0.1% w/v |
| Candida albicans | 11.2 * 10 ⁸ CFU/ml |
| Gardnerella vaginalis | 8.6 * 10 ⁸ CFU/ml |

Substance/Condition/ Micro-organism

Tested concentration

| | |
|-----------------------|--|
| Neisseria gonorrhoea | 10.6 * 10 ⁸ CFU/ml |
| Chlamydia trachomatis | 5.0 * 10 ³ TCID ₅₀ |
| HSV-1 | 5.0 * 10 ³ TCID ₅₀ |
| HSV-2 | 5.0 * 10 ³ TCID ₅₀ |

TRACEABILITY

The Actim PROM test is traceable to Behrings Institut's PP12 (IGFBP-1) preparation 307/323. This IGFBP-1 was purified from soluble extracts of human term placenta or decidua⁹ to prepare reference standards for a quantitative in-house method, IGFBP-1 IEMA, that is used in calibration.

DISPOSAL

Dispose of the contents according to the national and local law.

All patient specimens and used swabs, dipsticks and extraction buffer tubes should be handled and disposed of as potentially infectious material.

Materials used in Actim PROM kit and its components:

Paper: instructions for use.

Cardboard: Kit box.

Plastic: test pack pouch, extraction buffer tube.

Several (not to be recycled): unused dipstick, unopened extraction buffer tube, swabs, dipstick pouch.

When used in accordance with Good Laboratory Practice, good occupational hygiene and instructions for use, the reagents should not present a hazard to health.

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1. Silver, Thomas Lee and Helayne. Etiology and Epidemiology of Preterm Premature Rupture of the Membranes. 2001. Clinics in Perinatology. Vol. 28, Issue 4.
2. Rutanen E-M, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clinica Chimica Acta (1993) 214: 73–81.
3. Rutanen E-M et al. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clin Chim Acta (1996) 253:91-101.
4. Ragosch V et al. Insulin like growth factor binding protein 1 (IGFBP-1) and fetales Fibronectin in der Diagnostik eines vorzeitigen Blasensprunges. GebFra (1996) 56:1-6.
5. Jain K and Morris P. A clinical study to evaluate the usefulness of the MAST test in diagnosing pre-labour rupture of membranes. J Obstet Gynaecol (1998) 18:33-36.
6. Kubota T and Takeuchi H. Evaluation of insulin-like growth factor binding protein-1 as a diagnostic tool for rupture of the membranes. J Obstet Gynecol Res (1998) 24:411-417.
7. Erdemoglu E and Mungan T. Significance of detecting insulin-like growth factor binding protein-1 in cervicovaginal secretions: Comparison with nitrazine test and amniotic fluid volume assessment. Acta Obstet Gynecol Scand (2004) 83:622-626.
8. H.Bohn and W.Kreus. Isolierung und Charakterisierung eines neuen plazentaspezifischen Proteins (PP12). Arch.Gynecol (1980) 229: 279-291.

Explanation of symbols • Sýmbole selgitused • Vysvětlení symbolů • Symbolforklaring • Erklärung der Symbole • Ενεξήγηση των συμβόλων • Explicación de los símbolos • Symbolien selitykset • Explication des symboles • Szimbólumok magyarázata • Skýring á táknum • Spiegazione dei simboli • Simbolių aiškinimas • Simbolu skaidrojums • Uitleg van symbolen • Symbolforklaringer • Objaśnienie symboli • Explicação dos símbolos • Explicația simbolurilor • Описание символов • Vysvetlivky symbolov • Razlaga simbolov • Symbolförklaring • Sembollerin açıklamaları

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SWAB

Puritan Medical Products Company
31 School Street, Guilford, Maine 04443-0149, USA

2797

EC REP

EMERGO EUROPE Arnhem, The Netherlands

Actim[®] PROM

Actim Oy
Klovinpellontie 3, FI-02180 Espoo, Finland
Tel. +358 9 547 680
www.actimtest.com
actim@actimtest.com

© Actim Oy. All rights reserved.
Actim[®] is a trademark of Actim Oy.

AOACE30831-4 2025-07-30