

# Actim<sup>®</sup> Partus

A rapid diagnostic test  
to help assess the risk  
of preterm birth.



[actimtest.com](http://actimtest.com)

# Actim<sup>®</sup> Partus

## Technical specifications

Test principle	Qualitative immuno-chromatographic test
Target protein	phIGFBP-1 (phosphorylated form of insulin-like growth factor binding protein-1)
Specimen type	Cervical secretions
Specimen collection	Swab specimen from cervical os. during speculum examination
Reading time	5 minutes
Interpretation	Visual
Cut-off	10 µg/l in the extracted specimen
Test components	1 swab, 1 tube of extraction buffer, 1 dipstick
Number of tests/kit	10 tests (31931ETAC)
Storage	+2... +25°C



## Actim Partus Controls kit (31900ETAC) for quality control

Kit components	1 negative control (31980ETAC) 1 low positive control (31981ETAC) 1 high positive control (31982ETAC) 1 reconstitution solution (30889ETAC)
Storage	+2... +25°C

## Product information

Product	Product code
Actim Partus, 10 tests	31931ETAC
Actim Partus Controls	31900ETAC

**Actim Oy**  
– part of Medix Biochemica Group  
Klovinpellontie 3, FI-02180 Espoo, Finland  
actim@actimtest.com  
[www.actimtest.com](http://www.actimtest.com)

